

Convaid Trekker - HCPCS E1234
Canadian Order Form

Account #: Purchase Order #: Date:
Contact: Mark For:

Bill To: Ship To:
Address: Address:
City/State/Zip:
Phone: City/State/Zip:
Email: Phone:
Trekker Standard Features:
• TR14 Base Accepts TR12 & TR14 Seating
• Tilt-in-Space (-10° to +45°)
• 180° Reversible Seating Configuration
• 90° Fixed Front Rigging
• Single Foot Activated Wheel Lock Bar
• 7.5" x 2" Front and 11.5" x 2.5" Rear Solid Tires
• One Piece Height Adjustable Push Handle (Inside Mount)
• Silver Powder Coated Frame
• 90° Fixed Footplate
• Quick Release Rear Wheels
• Seat to Back Angle (80° to 110°)
• Depth Adjustable Seat Pan
• Made in the USA
Standard Models:
☐ Trekker TR14 - 904661

Trekker Seating Module
Back Module Options - Select one below:
☐ Back Canes Only (12)¹
☐ Back Canes Only (14)¹
☐ Height Adjustable Solid Back with Cushion (12) (Back Height Adjustment Range 13" - 17")¹
☐ Height Adjustable Solid Back with Cushion (14) (Back Height Adjustment Range 17" - 25")¹
Back Cushion:
☐ Back Cushion - Short (12: 12x18)
☐ Back Cushion - Short (14: 14x20)
☐ Back Cushion - Tall (14: 14x23)
Seating Module Options (See chart on last page of order form):
☐ TR12 Seat Depth Position 1
☐ TR12 Seat Depth Position 2
☐ TR12 Seat Depth Position 3
☐ TR14 Seat Depth Position 1
☐ TR14 Seat Depth Position 2
☐ TR14 Seat Depth Position 3
☐ TR14 Seat Depth Position 4
Upholstery Color:
☐ Panther Black
☐ Royal Blue
☐ Chocolate Brown
☐ Forest Green
☐ Sassy Purple
☐ Princess Pink
☐ Cherry Red
☐ Aqua Blue
☐ Apple Green
Extra Back Cushion Options\*
☐ 1" Open Cell Foam Insert Short Back (12)
☐ 1" Open Cell Foam Insert Short Back (14)
☐ 1" Open Cell Foam Insert Tall Back (14)
☐ 1/2" Soft Viscoelastic Foam Insert Short Back (12)
☐ 1/2" Soft Viscoelastic Foam Insert Short Back (14)
☐ 1/2" Soft Viscoelastic Foam Insert Tall Back (14)
☐ Lumbar Support
Contoured Seat Cushion\*\*
☐ Contoured Seat Cushion - Short (12: 12x11)
☐ Contoured Seat Cushion - Short (14: 14x13)
☐ Contoured Seat Cushion - Long (12: 12x14)
☐ Contoured Seat Cushion - Long (14: 14x16)
Planar Seat Cushion\*\*
☐ Planar Seat Cushion - Short (12: 12x12)
☐ Planar Seat Cushion - Short (14: 14x14)
☐ Planar Seat Cushion - Long (12: 12x15)
☐ Planar Seat Cushion - Long (14: 14x17)
Extra Seat Cushion Options\*
☐ 1" Open Cell Foam Insert Short Seat (Planar Cushion Required) (12)
☐ 1" Open Cell Foam Insert Short Seat (Planar Cushion Required) (14)
☐ 1" Open Cell Foam Insert Long Seat (Planar Cushion Required) (12)
☐ 1" Open Cell Foam Insert Long Seat (Planar Cushion Required) (14)
☐ 1/2" Soft Viscoelastic Foam Insert Short Seat (Planar Cushion Required) (12)
☐ 1/2" Soft Viscoelastic Foam Insert Short Seat (Planar Cushion Required) (14)
☐ 1/2" Soft Viscoelastic Foam Insert Long Seat (Planar Cushion Required) (12)
☐ 1/2" Soft Viscoelastic Foam Insert Long Seat (Planar Cushion Required) (14)

¹Manufacturer Installation Required.
²Fast tab is only added to Cushion if 3-Pt Belt is added. New Cushion is required if 3-Pt Belt is added after-market.
\*All seat cushions must have fast tab for aftermarket 3-Pt Belt and H-Harness.

Transit Required Items^
☐ Transit Option (Manufacturer Installation Required)
☐ Using other manufacturer's Head Support
^Requires Head Support and Pelvic Positioning Belt

Transit Options
☐ Q'Straint Transit Lap Belt
☐ Sure-Lok Transit Lap Belt

Table with 2 columns: Model Number, Description. Rows include TR14 models 904768 through 904719.

Table with 2 columns: Model Number, Description. Rows include TR14 models 904576 and 995681.

Table with 2 columns: Model Number, Description. Rows include TR14 models 903094 and 903232.

Pelvic Positioning

- ☐ 2-Point Pelvic Positioning Belt (Size 1)
- ☐ 3-Pt Positioning Belt (Trekker Seat Cushion Required)
- 4-Point Pelvic Positioning Belt by R82 (Hardware Required)**

☐ W: 9-10.25" X D: 21.5" (Size 0)

☐ W: 10.25-11.25" X D:21.5" (Size 1)

☐ W: 12-13.5" X D: 21.5" X H: 19.5" (Size 2)

☐ W: 14-15.5" X D: 21.5" X H: 19.5" (Size 3)

☐ R82 Hardware Kit (Required for R82 Belt)
- Left Width, Height Adj. Lateral Pelvic Support (Hip Guide w/out Armrest)**

☐ Small - 4.5" - 5.5" (H)\*

☐ Large - 5.25" - 6.25" (H)\*
- Right Width, Height Adj. Lateral Pelvic Support (Hip Guide w/out Armrest)**

☐ Small - 4.5" - 5.5" (H)\*

☐ Large - 5.25" - 6.25" (H)\*
- Left Width, Adj Pelvic Support w/Height Adj. Armrest (Hip Guide w/ Armrest) (n/a w/Tray)**

☐ Small - 5"-6"(H) (12)\*

☐ Small - 5"-6"(H) (14)\*

☐ Large - 7"-8" (H)\*
- Right Width, Adj Pelvic Support w/Height Adj. Armrest (Hip Guide w/ Armrest) (n/a w/Tray)**

☐ Small - 5"-6"(H) (12)\*

☐ Small - 5"-6"(H) (14)\*

☐ Large - 7"-8" (H)\*

\*Height is reduced by 1" when used with Planar Seat Cushion or 2" if used with Contoured Seat Cushion.

Foot & Leg Positioning

- ☐ Angle Adjustable Elevating Leg Rest (90° - 180°) Pair (12)<sup>1</sup>
- ☐ Angle Adjustable Elevating Leg Rest (90° - 180°) Pair (14)<sup>1</sup>
- ☐ Individual Angle Adjustable Footplates - Pair (Elevating Leg Rest Required) (12)<sup>1</sup>
- ☐ Individual Angle Adjustable Footplates - Pair (Elevating Leg Rest Required) (14)<sup>1</sup>
- ☐ One Piece Angle Adjustable Flip Up Footplate (-15° - +15°) (12)
- ☐ One Piece Angle Adjustable Flip Up Footplate (-15° - +15°) (14)
- ☐ Full Padded Footbox (One-Piece Footplate Required) (12)
- ☐ Full Padded Footbox (One-Piece Footplate Required) (14)
- Medial Knee Abductor**

☐ Small

☐ Medium
- Foot Positioners (Pair)\*\***

☐ Small 7" - 10"

☐ Medium 8" - 11"

☐ Large 9" - 12"
- Therafin Shoe Holders w/Padded Straps (Pair)\*\***

☐ Small

☐ Medium
- ☐ Calf Panel (E.I.) H: 3.5" ^^
- ☐ Calf Panel (Standard) H: 6"
- \*\*Cannot be used during transit.**
- <sup>1</sup>Manufacturer Installation Required.**
- ^^Use with Seat-to-Footplate range 4" - 7".**

TR14
903931
904017T
905148
905149
905150
905151
905309
904935
904937
904936
904938
904611
904613
904927
904612
904614
904928

TR14
905046
905047
905049
905050
904568
904596
904955
904956
904557
904591
903415
903416
903417
904679
904680
905370
904580

Trunk Positioning

- ☐ H-Harness w/ Padded Covers^^
- Neoprene Vest by R82 (Hardware Required)**

☐ W: 7" x D: 5.5" x H: 4" (Size 1)

☐ W: 6.5" x D: 9" x H: 6.75" (Size 2)

☐ W: 11.25" x D: 7.75" x H: 7" (Size 3)

☐ R82 Hardware Kit (Required for R82 Vest)
- Chest Belt, Upholstered by R82 (Hardware Required)**

☐ W: 11.25" x D:7.5" x H: 9.25-12" (Size 1)

☐ W: 14" X D: 10.25" X H: 10.25-12.75" (Size 2)

☐ W: 16.25" X D: 12.5" X H: 10.25-12.75" (Size 3)

☐ W: 20.25" X D: 16" X H: 11.25-14.5" (Size 4)

☐ R82 Hardware Kit (Required for R82 Belt)
- Cross Vest by R82 (Hardware Required)**

☐ W: 4.5" X D: .5" X H: 11" (Size 0)

☐ W: 5" X D: .5" X H: 12" (Size 1)

☐ W: 5.2" X D: .5" X H: 12.75" (Size 2)

☐ W: 6.25" X D: 4.75" X H: 15.25" (Size 3)

☐ W: 7.5" X D: 6.75" X H: 18.25" (Size 4)

☐ R82 Hardware Kit (Required for R82 Vest)
- Left - Width, Height Adj. Lateral Thoracic Support - Fixed Laterals**

☐ Size 1: 4" (D) x 4.5" (H)\*\*\*\*^

☐ Size 2: 5" (D) x 4.5" (H)\*\*\*\*^

☐ Size 3: 6" (D) x 4.5" (H)\*\*\*\*^
- Right - Width, Height Adj. Lateral Thoracic Support - Fixed Laterals**

☐ Size 1: 4" (D) x 4.5" (H)\*\*\*\*^

☐ Size 2: 5" (D) x 4.5" (H)\*\*\*\*^

☐ Size 3: 6" (D) x 4.5" (H)\*\*\*\*^
- Left - Height Adj. Lateral Thoracic Support - Swing Away Laterals**

☐ Size 1: 4" (D) x 2.75" (H)\*

☐ Size 2: 5.75" (D) x 4.25" (H)\*

☐ Swing Away Adj. Hardware (Required for S/A Laterals)
- Right - Height Adj. Lateral Thoracic Support - Swing Away Laterals**

☐ Size 1: 4" (D) x 2.75" (H)\*

☐ Size 2: 5.75" (D) x 4.25" (H)\*

☐ Swing Away Adj. Hardware (Required for S/A Laterals)
- \*Pad Size.**
- ^^Requires the Height Adjustable Solid Back System and 3-Pt Belt.**
- \*\*\*\*Requires the Height Adjustable Solid Back System.**
- ^Depth of Laterals is reduced by 1" with back upholstery.**

TR14
904695
905152
905153
905154
905309
905155
905156
905157
905158
905309
905176
905177
905178
905179
905180
905309
904760
904876
904877
904761
904886
904887
904864
904866
904884
904865
904867
904885

## TR14

- |        |
|--------|
| 904583 |
| 904607 |
| 904584 |
| 904608 |
|        |
| 904950 |
| 904951 |
| 904952 |
|        |
| 904948 |
| 904947 |
|        |
| 905282 |
| 905173 |
| 905288 |
| 905174 |
|        |
| 905284 |
| 905286 |
| 905175 |
| 905290 |
|        |
| 905281 |
| 905170 |
| 905289 |
| 905291 |
|        |
| 905285 |
| 905287 |
| 905172 |
| 905171 |
| 904946 |

905053
904688

- TR14**

- |        |
|--------|
| 905366 |
| 905367 |
| 905368 |
| 905369 |
| 905363 |

## TR14

- |        |
|--------|
| 904117 |
| 904118 |

## TR14

- |        |
|--------|
| 904566 |
| 904565 |
| 904671 |
| 904595 |

**TR14**

904572
904599
904690
904691
904573
904600
904556
904590
905361
904586
904610
905310
904689
904649
904605
904606
904531
904790
904791
904788
904789

- |        |
|--------|
| 904572 |
| 904599 |
| 904690 |
| 904691 |
| 904573 |
| 904600 |
| 904556 |
| 904590 |
| 905361 |
| 904586 |
| 904610 |
| 905310 |
| 904689 |
| 904649 |
| 904605 |
| 904606 |
| 904531 |
| 904790 |
| 904791 |
| 904788 |
| 904789 |

\*\*\*The Trekker Transit option includes four manufacturer installed red wheelchair transportation anchors. The Trekker Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are TR12T 75lbs/34kg and TR14T 110lbs/50kg.

The Convoid Trekker Mobility Base was tested with the Trekker Seating Module equipped with a head support and wheelchair anchored pelvic lap belt. Do not use the foot positioners on the Convoid Trekker during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is or can being made.

Furthermore, when being transported in a vehicle it is Convoid's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convoid cannot recommend a transit system.

For more information on RESNA WC-19 and WC-20 standards please visit <http://www.rercwts.pitt.edu/WC19.html>.

*HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.*

**Standard Shipping Rates: Ground \$55 • Expedited \$114 (Expedited shipping includes brokerage service and fee)**  
2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166  
Prices in CAD & subject to change, ex works Torrance, California.


www.convoid.com

Effective: 12/6/2018

Rev 9

**Seat Depth Adjustable Positions:**

Seat Depth Growth				
Position:	1	2	3	4
TR12:	7.5" to 9.5"	9.5" to 11.5"	11.5" to 13.5"	N/A
TR14:	9.5" to 11.5"	11.5" to 13.5"	13.5" to 15.5"	14.5" to 16.5"

**Seating Module Specifications:****Convaid Trekker Seating Module Specifications**

Seat width with accessories	TR 12 (Size 1)	TR 14 (Size 2)
Trunk width between fixed laterals	8.5-12"	10.5-14"
Trunk width between swing-away laterals	6.5-12"	8.5-14"
Hip width between hip guides without armrest	6-11.5"	6.5-14"
Hip width between hip guides with armrest	6.25-11.5"	6.5-14"
Seat to footplate dimensions	TR 12 (Size 1)	TR 14 (Size 2)
Fixed 90° with one piece footplate	5.5-11.5"	5.5-11.5"
Fixed 90° with one piece footplate (flipped)	3.25-9.25"	3.25-9.25"
ELR with one piece footplate	4.75-10.75"	6.25-12.25"
ELR with one piece footplate (flipped)	2-8"	4.25-10.25"
ELR with individual footplates	5-11"	7-13"
ELR with individual footplates (flipped)	2-5"	4-7"