



Convaid Trekker[®]

Our newest generation of Trekker wheelchair: Lightweight, compact-folding tilt-in-space wheelchair.



 **Convaid** | **R82** 
a part of Etac

Junior | Seating


CORDURA CUSHION COLORS



*Positioning options & accessories sold separately



Trekker Accessories



1 Extended Headrest Canopy
With or Without Windows

2 Headrest Cover (Canopy)

3 Mosquito Net
(Canopy Required)

4 Rain Cover
(Canopy Required)

5 Clear Upper Extremity
Support Tray

6 Padded Tray Cover
(Clear Tray Required)

7 IV Pole

8 Oxygen Tank Holder

9 Ventilator/Suction
Machine Hard Tray

10 Medical Necessity Basket

11 Under Seat Storage Basket

12 LTV Ventilator Holder

13 Caster Locks

14 Curb Tipper

15 Rear Anti-tip Tubes


















16 Transit Option



- | | | | |
|---|--|---|---|
| 17 Anatomic Headrest by R82 | 18 Adjustable Wing Head Support | 19 Height Adjustable Curved Head Support | 20 Utility Bag |
| 21 Chest Belt by R82 | 22 Cross Vest by R82 | 23 Neoprene Vest by R82 | 24 H-Harness with Padded Covers |
| 25 Width, Height, Adj. Lateral Thoracic Support - Fixed Laterals | 26 Height Adj. Lateral Thoracic Support - Swing-Away Laterals | 27 Width Adjustable Lateral Pelvic Support With or Without Armrest | 28 Medial Knee Abductor (Pommel) |
| 29 Contoured Back Cushion | 30 Contoured Seat Cushion
• 1" Open Cell Foam Insert (option)
• 1/2" Soft Viscoelastic Foam Insert (option) | 31 Planar Seat Cushion
• 1" Open Cell Foam Insert (option)
• 1/2" Soft Viscoelastic Foam Insert (option) | 32 A. 2-Point Positioning Belt
B. 3-Point Positioning Belt
C. 4-Point Positioning Belt |
| 33 Angle Adjustable Elevating Leg Rest with one piece footplate | 34 Individual Angle Adjustable. Footplates (ELR Required) | 35 Shoe Holders With Padded Straps (Therapin)
Footpositioners (not shown) | 36 Footbox
Calf Panel (not shown) |











Convoid Trekker Seat

Specifications

Seat Size	 1 (12")	 2 (14")
 Seat width	11.75"	13.75"
 Seat depth range	7.5" - 11.5"	9.5" - 16.5"
 Back height	13" - 17"	17" - 26"
 Hip width between hip guides	6.5"-11.5"	6"-14"
 Trunk width between fixed laterals	8.5"-12"	10.5"-14"
 Trunk width between swing-away laterals	6.5"-12"	8.5"-14"
 Seat to Footplate (one-piece footplate)	4" - 12" 3" - 11"*	4" - 12" 3" - 11"*
 Seat to Footplate (individual footplates) ^Δ	4.75" - 10.5" 1" - 5"*	7" - 13" 1" - 6.25"*
 Back angle adjustment	80° - 110°	80° - 110°
 Optional recline adjustment	80° - 170°	80° - 170°
 Tilt adjustment	-5° to 45°	-5° to 45°
 Headrest extension	8"	8"
 Shoulder strap height	14.5" - 20.5"	15" - 23"
 Footplate adjustment angle	±15°	±15°
 Weight of seating system ^Δ	10 lbs	11.5 lbs

Convoid Trekker Frame

Specifications

Frame Size	 1 (14")
 Folded length	39"
 Folded width	26"
 Folded height	13.5"
 Push handle adjustments	29" - 43"
 Weight of frame	25 lbs
 Overall height	37.5"
 Overall length	42.5"
 Overall width	26"
 Weight capacity regular/transit	110/110 lbs

*Reduced height when footplates are flipped
^ΔWeight of chair is without front rigging