

# Customized upper body support improves independence



## PivotFit™ Shoulder Harness

Promotes upright posture for improved breathing, head control and stability through appropriate shoulder positioning

- Comfortable support from body-contouring shape
- Individualized fit from high-strength swivel buckle
- Available in Dynamic or Non-Stretch styles



**Dynamic**  
Encourages movement to reduce tone



**Non-Stretch**  
Reduces movement to address positioning needs



## Stayflex™ Chest Support

Actively stabilizes the trunk while encouraging shoulder movement for long-term stability and performance

- Patented dual-zone elastic construction virtually eliminates upward shifting of chest support when the user leans forward
- Unique materials provide an optimal balance between stretch and resistance
- Standard and narrow styles available with and without zipper



**Standard**  
Comfortable with a controlled amount of stretch



**Narrow**  
Comfortably fits women and users with smaller frames



## Dynamic Monoflex™ Chest Support

Promotes free arm movement with minimal yet strong abdominal/chest support

- High-strength elastic padding comfortably stabilizes the trunk without over-stretching
- High-strength swivel buckles with double-tapered padding distribute pressure across the chest
- Available in Underarm Release with front-pull adjustment (above) or Center Release (below) styles



**Monoflex™ Center Release**  
Reduces bulk under the arms and fits under lateral trunk pads for added stability



## Trimline Shoulder Harness

Provides aggressive shoulder retraction with maximum freedom of movement across the chest

- Fits women comfortably and modestly with tapered pads and a buckle-free front
- Independent pads for versatile configurations: Crossed at the sternum (below) or backpack style (above)



\*Always use an anterior positioning device with a properly fitted pelvic support belt.

## 1 SELECT PULL-STYLE



**Rear-Pull**  
Tightens by pulling down and behind the user. Ideal for thicker backrests



**Front-Pull**  
Tightens by pulling forward, in front of the user



**Front-Pull Adjustable length**  
Top straps permit mounting to any type of backrest or mounting bracket

## 2 SELECT END-FITTING



**-B1 Flat-Mount™**  
Versatile, lowest profile

OR



**-B2 Cinch-Mount™**  
Time saving, precise, easy adjustment



**FS032-2 Cam Buckle, 1/4" Hole**  
For frequent length adjustment on Stayflex, PivotFit and Trimline Harnesses. Mount to solid backrest. Sold separately.