

Aluminum Crutches with Accessories



- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable accessories are latex free
- Weight Capacity 350 lbs on aluminum crutch (excluding pediatric)

PRODUCTS	Description
10400-8	8 pr/cs
10401-8	8 pr/cs
10402-8	8 pr/cs
10416-1	1 pr/cs

SPECIFICATIONS	Description
Weight Limit (10400 10401 10402)	350 lbs.
Weight Limit (10416)	175 lbs.
Warranty	Limited Lifetime