

Scooter Accessories



Lockable Front Box

Mounts onto the front of the tiller to secure valuables. Features a 10 lb./4.5 kg capacity.

ACC260 (silver), ACC270 (blue), 1551638 (red),

1537710 (slate grey)

Leo, Pegasus, Comet & Comet HD

1637910 (charcoal grey)

Pegasus^{METRO}, Pegasus^{PRO}, Comet^{PRO} & Comet^{ULTRA}



Lockable Rear Storage Box

Mounts into the rear accessory ports of the scooter for extra storage space.

1637913

Pegasus^{METRO}, Pegasus^{PRO}, Comet^{PRO} & Comet^{ULTRA}



Extra Large Front Basket

Mounts onto the front of the tiller to provide three times the basket space versus the standard composite basket.

10.5" H x 17.75" W x 10.5" D

1637915

Pegasus^{METRO}, Pegasus^{PRO}, Comet^{PRO} & Comet^{ULTRA}



Rear Basket

Large basket is ideal for storing extra items. 18" W x 13" D x 7" H.

ACC230

Suitable for all scooters.

Colibri requires 1148158, scooter accessory bracket.

Comet HD/Comet^{ULTRA} require 1541502, scooter extension tube.



Front and Rear Steel Bumpers

Mounts securely onto the scooter, providing additional protection for the frame and shroud.

1528275 - Front - 3W Pegasus, Pegasus^{METRO}

1638965 - Front - 4W Pegasus^{METRO}, Pegasus^{PRO}

1528298 - Front - Comet, Comet^{PRO}, Comet^{ULTRA}

1638966 - Rear - Pegasus^{METRO}, Pegasus^{PRO}

1638967 - Rear - Comet^{PRO}, Comet^{ULTRA}



Deluxe Mirror Set

Stylized front mirrors for mounting on the left and right side of the control panel. Sold as a pair.

1637911

Suitable for all scooters except the Colibri.



Standard Mirror

Removable mirror is adjustable and can be mounted on the left or right side of the control panel. Sold individually.

ACC250

Suitable for all scooters except the Colibri.



Beverage Holders

RAM Level Cup™ clamp-on drink holder.

TAGHCUP (left) Suitable for all scooters.

Foldable Cup Holder. Installs onto tiller pocket.

1521594 (right)

Suitable for all scooters except the Colibri and Leo.