

Bariatric Steel Transport Chair



- Supports individuals up to 450 lbs
- Reinforced steel frame provides added support
- Dual, reinforced steel cross braces
- Available in blue or red frames
- 12" rear flat free wheels
- Comes with swing-away footrests
- Padded armrests provide additional comfort
- Heavy-duty, nylon reinforced upholstery with a back carry pouch
- Removable, reversible desk-length arms
- Companion activated wheel locks

PRODUCTS	Description
BTR20-B	Blue, 1/cs
BTR20-R	Red, 1/cs
BTR22-B	Blue, 1/cs
BTR22-R	Red, 1/cs

SPECIFICATIONS	Description
Overall Width (BTR20)	24.5" (open) 9" (closed)
Overall Width (BTR22)	26.5" (open) 11.5" (closed)
Overall Height	37"
Seat-To-Floor Height	19.5"
Back Height	16"
Seat Dimensions (BTR20)	20" (W) x 18" (D)
Seat Dimensions (BTR22)	22" (W) x 18" (D)
Overall Length w/Riggings	38"
Weight (BTR20 BTR22)	45 lbs. 49 lbs.
Carton Shipping Weight (BTR20 BTR22)	50 lbs. 54 lbs.
Weight Limit	450 lbs.
Warranty:	Limited Lifetime

PARTS	Description
STDS2K3003	1/pr
STDS2KBTR20B	1/ea
STDS2KBTR20S	1/ea
STDS00I07I	
STDSDSKVB	1/ea
STDS2K303R	1/ea
STDS2K303L	1/ea

PARTS	Description
STDS2K3015	
STDS2K303I	
STDS2K3002	1/bag
STDS2K3032	1/pr
STDS2K3000	1/ea
STDS2K3045	1/ea
STDS2K30II	
STDS2K3042	1/ea
STDSI0II	
STDS2K3009	1/ea
STDS3J240I	
STDS8M4651	1/ea
STDS2K3010	
STDSI026	
STDS2K305R	1/cv
STDS2K305B	1/cv
STDS1072-8	1/ea
STDS2K30I2	
STDS001071	1/ea
STDS1011	1/ea
STDS1026	1/ea
STDS2K3006F	
STDS2K3010	
STDS2K3011	
STDS2K3012	
STDS2K3015	
STDS2K3031	
STDS3J2401	1/ea
STDS3J4503	1/ea
STDS4S3006F	
STDS4S30LB	
STDS4S30LR	
STDS4S30RB	
STDS4S30RR	
STDS8M4601	1/ea
STDS8M4603	1/ea