



# Raz-CAT

## Compact Attendant Tilt

Rehab Shower Commode

The **Raz-CAT** (Compact Attendant Tilt) has all of the features of our Raz-AT but in a more compact size. For pediatrics, adolescents or smaller adults, the Raz-CAT provides a fit that is better suited to physically-smaller clients. To accommodate client growth, the Raz-CAT has adjustable seat depths and is available with 3 back-width options.



## FEATURES

Extensive range of options available

Flip-up padded arms

Unique oblong seat aperture accommodates seat depth adjustment

Exclusive to Raz, Ischial Pelvic Align System (IPAS), provides 2" (5 cm) of fore/aft or 1" (2.5 cm) of lateral seat adjustment

Stainless steel frame with 4" (10 cm) of height adjustment

Dual-locking (swivel and roll) casters



Tension-adjustable back upholstery

Available back frames:  
14", 16", 18" (36 cm, 41 cm, 46 cm) wide  
15", 18" (38 cm, 46 cm) tall

Backrest assembly adjusts fore/aft for seat depth change

Adjustable back angle

Locking gas struts to facilitate tilting

### ALSO AVAILABLE

- Custom seats
- Shown with optional laterals, hip adductors and heel loops

SPECIFICATIONS	Raz-CAT (Compact Attendant Tilt)
WEIGHT CAPACITY	250 lbs (115 kg)
FRAME WIDTH (seat and Back frame)	18" (46 cm)
BACK FRAME WIDTHS	14" (35 cm), 16" (40 cm), or 18" (46 cm)
OVERALL WIDTH (with standard armpad)	21.75" (55 cm), 23" (58 cm) with 18" (46 cm) wide Back frame
SEAT WIDTH (commode seat)	18" (46 cm)
BACK HEIGHTS (top of seat to top of upholstery)	15" or 18" (38 cm or 46 cm)
RECLINING BACKS (top of seat to top of upholstery)	15" or 18" (38 cm or 46 cm)
SEAT DEPTH (front of back post to front of seat)	11" - 16" (28 cm - 40 cm)
OVERALL LENGTH WHEN UPRIGHT	32" (81 cm)
CLEARANCE HEIGHT UNDER SEAT	16.75", 17.75", 18.75", 19.75", 20.75" (42.5, 45, 47.5, 50, 53 cm) with pan brackets 17.5", 18.5", 19.5", 20.5", 21.5" (44.5, 47, 49.5, 52, 54.5 cm) without pan brackets
ARMREST HEIGHT ABOVE SEAT (20" Adjusta-Back)	15" (38 cm) High Back: 5.5" - 8" in ½" increments (14 - 20 cm) 18" (46 cm) High Back: 6.5" - 9" in ½" increments (16.5 - 23 cm)
SEAT-TO-FLOOR HEIGHT	21.5" - 25.5" in 1" increments (55 - 65 cm)
PRODUCT WEIGHT	55 lbs (25 kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)

*All measurements taken with standard molded seat and standard arm pad.*

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

*Specifications, and options are subject to change without prior notice.*

### MANUFACTURED BY:

**RAZ Design Inc.** 22 Howden Road, Toronto, Ontario Canada,

M1R 3C7 T: (877) 720-5678 E: [info@razdesigninc.com](mailto:info@razdesigninc.com)

RAZSALG004r01 Rev. Date: 2017-02-27

Raz Design products are distributed in Australia by:

### Paragon Mobility

3/10 Tapnor Crescent, Brendale, QLD, 4500 Australia

Tel: 13000 652 382

Fax: (07) 3112 6125

[www.paragonmobility.com.au](http://www.paragonmobility.com.au)

