

## All Terrain Canes

HCPCS:



FIGURE A




FIGURE B




FIGURE C




FIGURE D



FIGURE E




FIGURE F

- Walk safely and confidently on virtually any surface and under most weather conditions (Figure A - B)
- The retractable tip has four stainless steel spikes to provide a sturdy grip on slippery surfaces
- No other installation required and provides the best traction and stability
- Other manufacturers' items require separate installation and are less safe
- To prevent damage to floors, gripping claw should only be used outdoors
- Available in black or bronze
- Height adjustable from 30" to 39"
- Retail packaged

SPECIFICATIONS	Description
Height	30" - 39"
Warranty Limited	Lifetime
Weight Capacity	250 lbs