

Tranquille

Glider

SENIOR LIVING
BRODA®




Gentle gliding motion



*Stimulates active patients and
reduces the need to wander*



*Auto-locking feature increases
safety and prevents falls*



*Comfort Tension Seating®
increases support and sitting
tolerance*

Patient, Facility, & Caregiver Benefits

- Provides necessary stimulation for memory care patients
- Increases resident safety and contentment
- Promotes proper spinal posture and pelvic positioning

UPHOLSTERY OPTIONS



SPECIFICATIONS

Weight Capacity	250 lbs.* (114 kg)
Wheels	Rear transport wheel, wheels for non-patient transport**
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Common Accessories	Padded armrests, full or half tray, adjustable armrests, contoured headrest
Standard Padding	Seat and back pad vinyl cover

*Adjustments available to accommodate larger patients.
**The Tranquille Auto-Locking Glider is not a transport chair.

Model	Seat width	Seat depth	Seat height	Back height
Small low-back 100-10AL	18"	17"	16"	26"
Low-back 100-15AL	20"	17"	18"	26"
High-back 100-20AL	20"	17"	18"	32"

FRAME OPTIONS

