



Toilet Safety Frame With Height and Width Adjustable Arms



FIGURE A



FIGURE B

- Anodized aluminum is sturdy and lightweight
- Plastic armrests provide additional comfort and support
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions) (Figure A)
- Powder-coated aluminum bracket easily attaches frame to toilet
- Bracket secures safety frame in place (Figure B)

PRODUCTS

Description

12001-4	Assembled, 4/cs
12001KD-1	Knocked Down, 1/cs

SPECIFICATIONS

Description

Width Between Arms (12001-4)	18.5" (smallest) 22.75" (largest)
Width Between Arms (12001KD-1)	18" (smallest) 22" (largest)
Adjustable Arm Height	25.5" (lowest) 29.5" (highest)
Carton Shipping Weight (12001-4)	18 lbs.
Carton Shipping Weight (12001KD-1)	5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

12001BR-1	Toilet Safety Frame Bracket, 1/ea
-----------	-----------------------------------