



Pegasus^{METRO} 3-Wheel & 4-Wheel Outdoor

Model Numbers:

	3-Wheel	4-Wheel
Blue	PEGASUS3M-B	PEGASUS4M-B
Red	PEGASUS3M-R	PEGASUS4M-R
Silver	PEGASUS3M-S	PEGASUS4M-S
White	PEGASUS3M-W	PEGASUS4M-W
Green	PEGASUS3M-G	PEGASUS4M-G

The Pegasus^{METRO} three and four-wheel scooters are equipped with numerous features designed to make your outdoor driving experience an enjoyable one. Proprietary Speed Reduction Technology ensures safe maneuvering around corners, and an intuitive control panel with visual and audible signals combine with programmable electronics to make the Pegasus^{METRO} exceptionally easy to drive.

Features

- Invacare's Speed Reduction Technology (SRT).
- Ergonomic design to encourage better alignment of the head, shoulders and elbows for a more relaxed driving position.
- Infinite adjustment of the tiller between 70° and 90°.
- 4-pole motor provides reliable performance and a 12 km/h top speed.
- Full front and rear lighting system with automatic brake lights, turning signals and super bright, dual bulb headlights.
- Intuitive control panel with domed buttons, LED indicators and programmable audible signals for concise, easy operation.
- Contoured captain's chair is height adjustable and offers enhanced multi-position recline, fore/aft adjustment and armrest height, width and angle adjustment capabilities.
- Programmable electronics.
- Reinforced front and rear bumpers for improved durability.
- Adjustable rear and front spring suspension provide smooth travel on uneven terrain.
- 11" black tires with polished aluminum rims are capable of climbing obstacles of up to 5".
- Two-step brake release lever prevents accidental disengagement out of drive mode.
- Motor and electronics that are protected from corrosion and water, to safeguard the vital workings of the scooter.
- Dual rear accessory ports.
- 300 lbs. weight capacity.

Shroud Options



Classic Colour Shroud Options

Pegasus^{METRO}, Pegasus^{PRO}, Comet^{PRO} & Comet^{ULTRA}



Jasper Red



Jade Green



Marble White



Sapphire Blue



Diamond Silver

Premium Colour Shroud Options

Pegasus^{PRO} & Comet^{PRO}



Crystal Copper



Iconic Chrome



Electric Blue



Luxor Beige



Carmine Red

STYLE IS FOREVER. INVACARE SCOOTERS.

Scooter Accessories



Lockable Front Box

Mounts onto the front of the tiller to secure valuables. Features a 10 lb./4.5 kg capacity.

ACC260 (silver), ACC270 (blue), 1551638 (red),

1537710 (slate grey)

Leo, Pegasus, Comet & Comet HD

1637910 (charcoal grey)

Pegasus^{METRO}, Pegasus^{PRO}, Comet^{PRO} & Comet^{ULTRA}



Lockable Rear Storage Box

Mounts into the rear accessory ports of the scooter for extra storage space.

1637913

Pegasus^{METRO}, Pegasus^{PRO}, Comet^{PRO} & Comet^{ULTRA}



Extra Large Front Basket

Mounts onto the front of the tiller to provide three times the basket space versus the standard composite basket.

10.5" H x 17.75" W x 10.5" D

1637915

Pegasus^{METRO}, Pegasus^{PRO}, Comet^{PRO} & Comet^{ULTRA}



Rear Basket

Large basket is ideal for storing extra items. 18" W x 13" D x 7" H.

ACC230

Suitable for all scooters.

Colibri requires 1148158, scooter accessory bracket.

Comet HD/Comet^{ULTRA} require 1541502, scooter extension tube.



Front and Rear Steel Bumpers

Mounts securely onto the scooter, providing additional protection for the frame and shroud.

1528275 - Front - 3W Pegasus, Pegasus^{METRO}

1638965 - Front - 4W Pegasus^{METRO}, Pegasus^{PRO}

1528298 - Front - Comet, Comet^{PRO}, Comet^{ULTRA}

1638966 - Rear - Pegasus^{METRO}, Pegasus^{PRO}

1638967 - Rear - Comet^{PRO}, Comet^{ULTRA}



Deluxe Mirror Set

Stylized front mirrors for mounting on the left and right side of the control panel. Sold as a pair.

1637911

Suitable for all scooters except the Colibri.



Standard Mirror

Removable mirror is adjustable and can be mounted on the left or right side of the control panel. Sold individually.

ACC250

Suitable for all scooters except the Colibri.



Beverage Holders

RAM Level Cup™ clamp-on drink holder.

TAGHCUP (left) Suitable for all scooters.

Foldable Cup Holder. Installs onto tiller pocket.

1521594 (right)

Suitable for all scooters except the Colibri and Leo.

Scooter Accessories



Mobile Phone Holder

Sturdy and secure bracket which adjusts to fit and hold most smartphones.

1637912

Suitable for all scooters except the Colibri.



USB Charger

Plugs into the scooters' charger port to enable the charging of one mobile device.

1637917

Pegasus^{METRO}, Pegasus^{PRO}, Comet^{PRO} & Comet^{ULTRA}



Safety Flag

6' bright orange flag makes the scooter easy to see from a distance. Comes complete with a secure mounting bracket.

ACC140

Suitable for all scooters.

Colibri requires 1148158, scooter accessory bracket.

Comet HD/Comet^{ULTRA} require 1541502, scooter extension tube.



Forearm Crutch and Cane Holder

Keeps a cane or forearm crutches stable and within reach during transit.

ACC200

Suitable for all scooters.

Colibri requires 1148158, scooter accessory bracket.

Comet HD/Comet^{ULTRA} require 1541502, scooter extension tube.



Walker Holder

Keeps a folding walker or rollator secure during travel.

ACC210

Suitable for all scooters.

Colibri requires 1148158, scooter accessory bracket.

Comet HD/Comet^{ULTRA} require 1541502, scooter extension tube.



Oxygen Tank Holder

Secure storage of up to a D-size oxygen tank for transit.

ACC220

Suitable for all scooters.

Colibri requires 1148158, scooter accessory bracket.

Comet HD/Comet^{ULTRA} require 1541502, scooter extension tube.



Windshield

Installs into the mirror ports by using either the mirrors or the hardware included.

1543160

Suitable for all scooters except the Colibri.

Scooter Accessories



Backrest Bag

Easily installs onto the scooter seat for additional storage.

PRC571720 (left)/ 1505731 (right)

Suitable for all scooters, except the Colibri and Leo.



Seat Suspension

Installs into the seat post for a softer ride.

1577212 (L)/1577213 (M)/1577214 (T) - Colibri

1541628 (L)/1541629 (M) - Leo

1640428 (L)/1640430 (M)/1640448 (T) - Comet^{PRO} to 201 lbs.

1640429 (L)/1640447 (M)/1640449 (T) - Comet^{PRO} to 285 lbs.

1640422 (L)/1640424 (M)/1640426 (T) - Pegasus^{Metro/Pro} to 201 lbs.

1640423 (L)/1640425 (M)/1640427 (T) - Pegasus^{Metro/Pro} to 285 lbs.



Scooter Cover

Heavy duty black nylon cover protects the scooter from the elements.

PRC571840 - Colibri & Leo

1518240 - Pegasus, Pegasus^{Metro/Pro}, Comet, Comet^{PRO}

Seat Cover

Protects seat from dirt rain and dirt.

1521684 (Not shown) - Pegasus, Pegasus^{Metro/Pro}, Comet, Comet^{PRO}



Enhanced Visibility Safety Vest

Highly reflective safety vest. One size.

SP1521592



Phone/Wallet/Comet HD

Pneumatic Bulb Arm Pocket

Attaches scooter seat armrest

12511

Suitable for all scooters.



Accessory Bracket

Attaches to the seat frame of the Colibri to accommodate one rear accessory.

1148158 (left)

Extension Tube

Extends the length of the Comet HD/ Comet^{ULTRA} accessory tube to accommodate the use of rear-mounted accessories.

1541502 (right)

Product Specifications

Overall Dimensions	Colibri 3W 12Ah/ 4W 12Ah/4W 18Ah	Leo 3W/4W	Pegasus 3W/4W METRO	Pegasus 4W PRO	Comet 4W PRO	Comet 4W ULTRA
Base Width	19.25"	23.6"	26"	25.5"	26"	26"
Max. climbable obstacle height	1.75"/1.75"/3"	2.4"/2.5"	3"	4"	4"	4"
Length	40"	47.5"	49"/50"	52"	57"	57"
Height	34.1"	36"	49.5"	49.5"	49.5"	50.75"
Usable Foot Space Length	18"/12"/12"	22"/13"	23"/13"	13"	13.5"	13.5"
Seat Dimensions						
Seat-to-Floor	14.5" - 16.5"	16.5" - 19.5"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"	17.25" - 20.25"
Width	17"	19"	20"	20"	20"	26"
Depth	17.5"	16"	18.5"	18.5"	18.5"	22.5"
Back Height	15.5"	18.5"	27" - 29"	27" - 29"	27" - 29"	32" - 35"
Back Angle	N/A	N/A	90° - 130°	90° - 130°	90° - 130°	95° - 130°
Arm Dimensions						
Length	11.5"	14"	14"	14"	14"	14"
Seat-to-Top of Arms	8"	8"	7" - 10"	7" - 10"	7" - 10"	5" - 7"
Width Between Arms	16" - 19"	23" - 27"	20" - 22.5"	20" - 22.5"	20" - 22.5"	20" - 22.5"
Recommended Battery	2 x 12 V, 12 Ah*/	2 x 12 V, 36 Ah*	2 x 12 V, 50 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah	2 x 12 V, 75 Ah
	2 x 12 V, 18 Ah*		Group 22	Group 24	Group 24	Group 24
Off-Board Charging System	2 Ah/2 Ah/5Ah	5 Ah	8 Ah	8 Ah	8 Ah	8 Ah

*included with scooter

Product Specifications

	Colibri 3W 12Ah/ 4W 12Ah/4W 18Ah	Leo 3W/4W	Pegasus 3W/4W METRO	Pegasus 4W PRO	Comet 4W PRO	Comet 4W ULTRA Pro
Speed	8 km/h (5 mph)	8 km/h (5 mph)	12 km/h (7.5 mph)	15 km/h (9.5 mph)	15 km/h (9.5 mph)	10 km/h (6 mph)
Estimated Travel Range*	11 km (7 miles) /	36 km (22 miles)	42 km (26 miles)	55 km (34 miles)	55 km (34 miles)	58 km (36 miles)
	11 km (7 miles) /					
	16 km (10 miles)					
Weight Limit (lbs.)	300**	350**	300**	350**	390**	485**
Maximum Incline Capability	6°	6°/8°	6°	10°	10°	10°
Tires						
Dimensions	200 x 50 / 200 x 50 / 210 x 65	9" x 3.5"	11" x 3.5" 4.10/3.50-5	12" x 3.5" 3.50-5	13" x 4" 4.00-5	13" x 4" 4.00-5
Tire Pressure	N/A	50 psi	50 psi	50 psi	50 psi	50 psi
Weight (lbs.)	91.6 / 97.4 / 108 (inc. batteries)	174/184 (inc. batteries)	114/125	172	185	212
Warranty Two-year limited warranty against manufacturing defects (excluding upholstered materials, padded materials and tires/wheels). Six months on batteries purchased from Invacare®. Manufacturer's warranty does not cover normal wear.						

*Travel range is an estimate only based on ISO 7176-2 standards determined by testing using a 100 kg/220 lb. dummy. Travel range will vary depending on user weight, speed, battery type, battery condition and charge level, weather, terrain and condition of the drive wheels.

**Weight limit is the maximum total weight of the rider, accessories and additional items the rider may be transporting. Example: If the weight limit of the scooter is 350 lbs. and additional items equal 25 lbs., the maximum rider weight can be 325 lbs.