




## Permobil Modular Systems Retail Packages



Buy together and save!



April 2018

To our valued clients,

These modular surfaces provide the optimal experience for each client. With just one mattress, there are multiple configurations available by repositioning the ROHO or foam sections\*.

The **Permobil Advantage**: Exceptional Products, Service and Price!

For more information, consult our complete product guide, price list, or contact customer service.

\*Substitute the High Resiliency Foam to GelFlex™ or Visco Foam and only pay the difference.

## Table of Contents

Acute Series.....	3
Original Snap Series - 36" .....	4
Original Snap Series - 39" .....	5
Adapt 2 Series - 36" .....	6
Adapt 2 Series - 39" .....	7

## ROHO® Mattress Sections

We have expanded our line of ROHO® compatible mattresses to give you even more options. These non-powered, adjustable, advanced pressure relieving mattress sections can be inserted into our: **Acute, Adapt 2, and Original Snap Series**.

These mattress sections are flame and fluid resistant and can be used in a number of configurations - whatever best suits your needs!



## Recovery 5 Incontinent Covers and Coverlets

All of our mattresses and inserts come standard with a Recovery 5 anti-microbial base cover which provides incontinent protection for up to four hours, and aids in keeping skin cool and dry. It is removable, machine-washable, fluid-proof and fire-retardant. This material helps wick away moisture and increases air flow and breathability.

**NOTE: Top covers are sold separately.** Top covers are available in the new Recovery 5 Ultra material which offers 4-way stretch, is machine washable, and provides up to 6X more incontinent protection. To order, simply add the code -R5U to the mattress part numbers, and only pay the difference for the upgraded material.

You can add a coverlet to your mattress which has stretchy corners and fits over the mattress (and inserts) like a regular fitted sheet. Each insert can also have its own removable zippered cover.



## Zippered Enclosure Cover

The OSS and Adapt 2 mattresses both have the option of adding a zippered enclosure cover in either the Recovery 5 or Recovery 5 Ultra material. The removable top cover with its 360° zipper must be ordered separately.



## Mattress Specifications and Warranty

Mattress Weight Capacity: 350 lb

Standard Mattress Size: 36" x 80" and 39" x 80"


Limited Warranty: **Two year complete coverage on mattress and covers**



**5MM mattresses  
are proudly made  
in Canada**

## To Place Your Order


Permobil Canada is the exclusive distributor for 5 Minds Mobility's mattresses. To order, please contact us at sales.ca@permobil.com or 800-265-9830. To learn more, visit our websites at www.permobil.ca and www.5minds-mobility.com.




## Acute Series Mattress Systems - 36" x 80"

The Acute mattress is our modifiable foam surface for the client at high risk of experiencing skin breakdown and may require the addition of ROHO® air components. Weight Capacity: 350 lb

All of our mattresses and inserts come with a Recovery 5 anti-microbial base cover. Top covers are sold separately.



Part Number	Description
AC-SYS1-PKG	ACUTE SYSTEM 1 PKG. INCLUDES 1 ROHO® MATTRESS SECTION (19" x 34")
AC-SYS2-PKG	ACUTE SYSTEM 2 PKG. INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 34")
AC-SYS3-PKG	ACUTE SYSTEM 3 PKG. INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 34"), MOVABLE DIVIDER
AC-SYS4-PKG	ACUTE SYSTEM 4 PKG. INCLUDES 3 ROHO® MATTRESS SECTIONS (19" x 34")
AC-SYS5-PKG	ACUTE SYSTEM 5 PKG. INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 34"), IMMOVABLE DIVIDER
AC-SYS6-PKG	ACUTE SYSTEM 6 PKG. INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 34") (T-SHAPE)



Acute system 1	Acute system 2	Acute system 3	Acute system 4	Acute system 5	Acute system 6
FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER, MOVABLE DIVIDER INCLUDED	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER, IMMOVABLE DIVIDER	FOAM BASE WITH RECOVERY 5 BASE COVER, T-SHAPE OPENING
1 x ROHO® MATTRESS SECTION (19" x 34")	2 x ROHO® MATTRESS SECTIONS (19" x 34")	2 x ROHO® MATTRESS SECTIONS (19" x 34")	3 x ROHO® MATTRESS SECTIONS (19" x 34")	2 x ROHO® MATTRESS SECTIONS (19" x 34")	2 x ROHO® MATTRESS SECTIONS (19" x 34")

### Cover Options - sold separately & new lower prices!

Part Number	Description
AC-COVERLET9	ACUTE - 9" FITTED ELASTIC RECOVERY 5 TOP COVER 36" x 80"
AC-COVERLET9-R5U	ACUTE - 9" FITTED ELASTIC RECOVERY 5 ULTRA TOP COVER 36" x 80"
AC-RC1934	ACUTE - RECOVERY 5 TOP COVER FOR ROHO® MATTRESS SECTION (19" x 34")

NOTE: Substitute the High Resiliency Foam to GelFlex™ or Visco Foam and only pay the difference.

Custom sizes available, contact customer service for further details.


## Original Snap Series Mattress System - 36" x 80"

The Original Snap mattress is our modifiable foam surface for the client at high risk of experiencing skin breakdown and may require the addition of ROHO® air mattress sections. Snaps hold the interchangeable inserts in place. Weight Capacity: 350 lb

All of our mattresses and inserts come with a Recovery 5 anti-microbial base cover. Top covers are sold separately.



Part Number	Description
SS-SYS1-36	ORIGINAL SNAP SERIES SYSTEM INCLUDES 1 ROHO® MATTRESS SECTION (19" x 30"), 1 FOAM WITH FEMALE SNAPS (30" x 30"), 1 FOAM WITH MALE SNAPS (30" x 30")
SS-SYS2-36	ORIGINAL SNAP SERIES SYSTEM INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 30"), 1 FOAM WITH SNAPS (30" x 30"), 1 FOAM WITH SNAPS (8" x 30"), 1 FOAM WITH SNAPS (4" x 30")
SS-SYS3-36	ORIGINAL SNAP SERIES SYSTEM INCLUDES 3 ROHO® MATTRESS SECTIONS (19" x 30"), 1 FOAM WITH SNAPS (20" x 30"), 1 FOAM WITH SNAPS (4" x 30")
SS-SYS4-36	ORIGINAL SNAP SERIES SYSTEM INCLUDES 4 ROHO® MATTRESS SECTIONS (19" x 30"), 1 FOAM WITH SNAPS (4" x 30")
SS-SYS5-36	ORIGINAL SNAP SERIES SYSTEM INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 30"), 1 FOAM WITH FEMALE SNAPS (30" x 30"), 1 FOAM WITH SNAPS (12" x 30")



System 1 / SS-SYS1-36	System 2 / SS-SYS2-36	System 3 / SS-SYS3-36	System 4 / SS-SYS4-36	System 5 / SS-SYS5-36
FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER
1 x ROHO® MATTRESS SECTION (19" x 30")	2 x ROHO® MATTRESS SECTIONS (19" x 30")	3 x ROHO® MATTRESS SECTIONS (19" x 30")	4 x ROHO® MATTRESS SECTIONS (19" x 30")	2 x ROHO® MATTRESS SECTIONS (19" x 30")
1 x SS-3030FS-F	1 x SS-3030FS-F	1 x SS-2030FS	1 x SS-430FS	1 x SS-3030FS-F
1 x SS-3030FS-M	1 x SS-830FS	1 x SS-430FS		1 x SS-1230FS
	1 x SS-430FS			

### Cover Options - sold separately & new lower prices!

Part Number	Description
AD2-TC3680	REMOVABLE ZIPPERED RECOVERY 5 TOP COVER 36" x 80"
AD2-TC3680-R5U	REMOVABLE ZIPPERED RECOVERY 5 ULTRA TOP COVER 36" x 80"
AC-COVERLET9	9" FITTED ELASTIC RECOVERY 5 TOP COVER 36" x 80"
AC-COVERLET9-R5U	9" FITTED ELASTIC RECOVERY 5 ULTRA TOP COVER 36" x 80"
COV-MATTSECT1610	RECOVERY 5 TOP COVER FOR ROHO® MATTRESS SECTION (19" x 30")

NOTE: Substitute the High Resiliency Foam to GelFlex™ or Visco Foam and only pay the difference.

Custom sizes available, contact customer service for further details.










## Original Snap Series Mattress System - 39" x 80"

The Original Snap mattress is our modifiable foam surface for the client at high risk of experiencing skin breakdown and may require the addition of ROHO® air mattress sections. Snaps hold the interchangeable inserts in place. Weight Capacity: 350 lb

All of our mattresses and inserts come with a Recovery 5 anti-microbial base cover. Top covers are sold separately.



Part Number	Description
SS-SYS1-39	ORIGINAL SNAP SERIES SYSTEM INCLUDES 1 ROHO® MATTRESS SECTION (19" x 34"), 1 FOAM WITH FEMALE SNAPS (30" x 32"), 1 FOAM WITH MALE SNAPS (30" x 32")
SS-SYS2-39	ORIGINAL SNAP SERIES SYSTEM INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 34"), 1 FOAM WITH SNAPS (30" x 32"), 1 FOAM WITH SNAPS (8" x 32"), 1 FOAM WITH SNAPS (4" x 32")
SS-SYS3-39	ORIGINAL SNAP SERIES SYSTEM INCLUDES 3 ROHO® MATTRESS SECTIONS (19" x 34"), 1 FOAM WITH SNAPS (20" x 32"), 1 FOAM WITH SNAPS (4" x 32")
SS-SYS4-39	ORIGINAL SNAP SERIES SYSTEM INCLUDES 4 ROHO® MATTRESS SECTIONS (19" x 34"), 1 FOAM WITH SNAPS (4" x 32")
SS-SYS5-39	ORIGINAL SNAP SERIES SYSTEM INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 34"), 1 FOAM WITH FEMALE SNAPS (30" x 32"), 1 FOAM WITH SNAPS (12" x 32")



System 1 / SS-SYS1-39	System 2 / SS-SYS2-39	System 3 / SS-SYS3-39	System 4 / SS-SYS4-39	System 5 / SS-SYS5-39
FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER
1 x ROHO® MATTRESS SECTION (19" x 34")	2 x ROHO® MATTRESS SECTIONS (19" x 34")	3 x ROHO® MATTRESS SECTIONS (19" x 34")	4 x ROHO® MATTRESS SECTIONS (19" x 34")	2 x ROHO® MATTRESS SECTIONS (19" x 34")
1 x SS-3032FS-F	1 x SS-3032FS-F	1 x SS-2032FS	1 x SS-432FS	1 x SS-3032FS-F
1 x SS-3032FS-M	1 x SS-832FS	1 x SS-432FS		1 x SS-1232FS
	1 x SS-432FS			

### Cover Options - sold separately & new lower prices!

Part Number	Description
AD2-TC3980	REMOVABLE ZIPPERED RECOVERY 5 TOP COVER 39" x 80"
AD2-TC3980-R5U	REMOVABLE ZIPPERED RECOVERY 5 ULTRA TOP COVER 39" x 80"
COV-MATTSECT3400	RECOVERY 5 TOP COVER FOR ROHO® MATTRESS SECTION (19" x 34")

NOTE: Substitute the High Resiliency Foam to GelFlex™ or Visco Foam and only pay the difference.

Custom sizes available, contact customer service for further details.

## Adapt 2 Series Mattress System - 36" x 80"

The Adapt 2 mattress with SRT is an innovative way to reduce the effects of pressure, friction and shear within the body of the mattress. The transfer sides and 4" permanent GelFlex™ foot bolster hold the inserts in place without the need for snaps. The surface inserts can be configured based on the individual client's needs. The removable top cover with its 360° zipper must be ordered separately. Weight Capacity: 350 lb


**Note:** Adapt 2 surface and inserts include a Fire Barrier sock fabric to meet full-scale fire retardant standards (not pictured).

All of our mattresses and inserts come with a Recovery 5 anti-microbial base cover. Top covers are sold separately.

Our newest offering!



Part Number	Description
AD2-SYS1-36	ADAPT 2 SERIES SYSTEM INCLUDES 1 ROHO® MATTRESS SECTION (19" x 30"), 1 FOAM INSERT (30" x 30"), 1 FOAM INSERT (19" x 30"), 1 FOAM INSERT (8" x 30")
AD2-SYS2-36	ADAPT 2 SERIES SYSTEM INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 30"), 1 FOAM INSERT (30" x 30"), 1 FOAM INSERT (8" x 30")
AD2-SYS3-36	ADAPT 2 SERIES SYSTEM INCLUDES 3 ROHO® MATTRESS SECTIONS (19" x 30"), 1 FOAM INSERT (19" x 30")
AD2-SYS4-36	ADAPT 2 SERIES SYSTEM INCLUDES 4 ROHO® MATTRESS SECTIONS (19" x 30")
AD2-SYS5-36	ADAPT 2 SERIES SYSTEM INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 30"), 1 FOAM INSERT (30" x 30"), 1 FOAM INSERT (8" x 30")




System 1 / AD2-SYS1-36	System 2 / AD2-SYS2-36	System 3 / AD2-SYS3-36	System 4 / AD2-SYS4-36	System 5 / AD2-SYS5-36
FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER
1 x ROHO® MATTRESS SECTION (19" x 30")	2 x ROHO® MATTRESS SECTIONS (19" x 30")	3 x ROHO® MATTRESS SECTIONS (19" x 30")	4 x ROHO® MATTRESS SECTIONS (19" x 30")	2 x ROHO® MATTRESS SECTIONS (19" x 30")
1 x AD2-3030FS	1 x AD2-3030FS	1 x AD2-1930FS		1 x AD2-3030FS
1 x AD2-1930FS	1 x AD2-830FS			1 x AD2-830FS
1 x AD2-830FS				

### Cover Options - sold separately & new lower prices!

Part Number	Description
AD2-TC3680	REMOVABLE ZIPPERED RECOVERY 5 TOP COVER 36" x 80"
AD2-TC3680-R5U	REMOVABLE ZIPPERED RECOVERY 5 ULTRA TOP COVER 36" x 80"
AC-COVERLET9	9" FITTED ELASTIC RECOVERY 5 TOP COVER 36" x 80"
AC-COVERLET9-R5U	9" FITTED ELASTIC RECOVERY 5 ULTRA TOP COVER 36" x 80"
COV-MATTSECT1610	RECOVERY 5 TOP COVER FOR ROHO® MATTRESS SECTION (19" x 30")

NOTE: Substitute the High Resiliency Foam to GelFlex™ or Visco Foam and only pay the difference.

Custom sizes available, contact customer service for further details.




## Adapt 2 Series Mattress System - 39" x 80"

The Adapt 2 mattress with SRT is an innovative way to reduce the effects of pressure, friction and shear within the body of the mattress. The transfer sides and 4" permanent GelFlex™ foot bolster hold the inserts in place without the need for snaps. The surface inserts can be configured based on the individual client's needs. The removable top cover with its 360° zipper must be ordered separately. Weight Capacity: 350 lb






**Note:** Adapt 2 surface and inserts include a Fire Barrier sock fabric to meet full-scale fire retardant standards (not pictured).

All of our mattresses and inserts come with a Recovery 5 anti-microbial base cover. Top covers are sold separately.

Our newest  
offering!



Part Number	Description
AD2-SYS1-39	ADAPT 2 SERIES SYSTEM INCLUDES 1 ROHO® MATTRESS SECTION (19" x 34"), 1 FOAM INSERT (30" x 33"), 1 FOAM INSERT (19" x 33"), 1 FOAM INSERT (8" x 33")
AD2-SYS2-39	ADAPT 2 SERIES SYSTEM INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 34"), 1 FOAM INSERT (30" x 33"), 1 FOAM INSERT (8" x 33")
AD2-SYS3-39	ADAPT 2 SERIES SYSTEM INCLUDES 3 ROHO® MATTRESS SECTIONS (19" x 34"), 1 FOAM INSERT (19" x 33")
AD2-SYS4-39	ADAPT 2 SERIES SYSTEM INCLUDES 4 ROHO® MATTRESS SECTIONS (19" x 34")
AD2-SYS5-39	ADAPT 2 SERIES SYSTEM INCLUDES 2 ROHO® MATTRESS SECTIONS (19" x 34"), 1 FOAM INSERT (30" x 33"), 1 FOAM INSERT (8" x 33")



System 1 / AD2-SYS1-39	System 2 / AD2-SYS2-39	System 3 / AD2-SYS3-39	System 4 / AD2-SYS4-39	System 5 / AD2-SYS5-39
FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER	FOAM BASE WITH RECOVERY 5 BASE COVER
1 x ROHO® MATTRESS SECTION (19" x 34")	2 x ROHO® MATTRESS SECTIONS (19" x 34")	3 x ROHO® MATTRESS SECTIONS (19" x 34")	4 x ROHO® MATTRESS SECTIONS (19" x 34")	2 x ROHO® MATTRESS SECTIONS (19" x 34")
1 x AD2-3033FS	1 x AD2-3033FS	1 x AD2-1933FS		1 x AD2-3033FS
1 x AD2-1933FS	1 x AD2-833FS			1 x AD2-833FS
1 x AD2-833FS				

### Cover Options - sold separately & new lower prices!

Part Number	Description
AD2-TC3980	REMOVABLE ZIPPERED RECOVERY 5 TOP COVER 39" x 80"
AD2-TC3980-R5U	REMOVABLE ZIPPERED RECOVERY 5 ULTRA TOP COVER 39" x 80"
COV-MATTSECT3400	RECOVERY 5 COVER FOR ROHO® MATTRESS SECTION (19" x 34")

NOTE: Substitute the High Resiliency Foam to GelFlex™ or Visco Foam and only pay the difference.  
Custom sizes available, contact customer service for further details.