

Molift RgoSling Ambulating Vest

For gait training, independent or with assistance.



Secure and comfortable gait training

Molift RgoSling Ambulating Vest is developed for safe and secure rehabilitation such as gait training, standing exercise as well as balance exercise.

Molift RgoSling Ambulating Vest has a double buckled wide padded belt that distributes pressure evenly over a larger surface. The padded groin straps limit the vest from sliding up and out of position.

Adopted features from professional climbing

Molift RgoSling Ambulating Vest has adopted features from professional climbing equipment to provide the best possible comfort and safety for the user. The self-tightening suspension loops prevent users from slipping or falling down.

The front buckle system provides increased possibilities for the user to more independently perform gait training or toileting, maybe even without assistance.

Accessory



Groin Straps

The groin straps are placed in an outward angle which minimizes risk of pinching when more weight is put on the vest.

Size	Item no.
XXS	1721300
XS	1721310
S	1721320
M	1721330
L	1721340
XL	1721350

Accessories

Groin Straps for Ambulating vest
Order two for a complete pair.

Size	Item no.
XXS–XS	1721700
S–M	1721720
L–XL	1721740

Material
Polyester,
PE-cellfoam

Recommended user weight

XXS	12–17 kg (26–37 lbs)
XS	17–25 kg (37–55 lbs)
S	25–50 kg (55–110 lbs)
M	45–95 kg (99–210 lbs)
L	90–160 kg (198–353 lbs)
XL	160–240 kg (353–529 lbs)

Washing Instruction



Colour code

