

rifton Compass Chair



permobil

Distributed in Canada by Permobil www.permobil.ca
Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944



Compass Chair (R330)

Compass Chair features



Unique built-in support

The contoured sides of the seat, together with the inside surfaces of the armrests, provide subtle lateral support that can be particularly helpful for students with Autism Spectrum Disorder, giving security and sensory boundaries.



The Compass Chair comes in five sizes.



Stability feet attach to the back legs to accommodate your more active, tip-prone students.



Adjusting the rear legs shorter than the front legs gives the chair a slight tilt, which can help some students relax.

Components



Seat belt

The optional seat belt clips into the seat when needed for safety or postural support.



Stability feet

The stability feet enlarge the chair's base of support. They are adjustable like the standard chair legs.



Adjustable legs

The legs provide 5 cm of height adjustment to ensure optimal seat height and allow for growth. Both chairs in this picture are size 3.

You asked; we delivered.



The return of an old friend

For years you've been asking for a simple supportive chair with armrests – like the old Compass Chair. Studying your suggestions and collaborating with therapists and teachers around the world, our team worked to develop the perfect (we hope!) inclusion chair.

While we kept the name, this is more than the reintroduction of an old product. The new chair has height adjustment, improved ergonomics and better support. We think it's the perfect chair for learning. Tell us if we got it right.

WWW.RIFTON.COM