

COMPANY _____

SHIP TO _____

DATE **MM/DD/YYYY** _____

NAME _____ PHONE _____

CONTACT NAME _____

ADDRESS _____

EMAIL _____

CITY _____ PROV _____ POSTAL CODE _____

MODEL SELECTION

- ☐ Z360 **RAZ-AT600** Heavy Duty Rehab Shower Commode Chair with stainless steel frame and components, 25° tilt range, four locking 6" (150mm) casters and commode pan, and the specifically engineered, Weight Assist Spring System (WASS).

Model Code		Seat Type	Frame Width	Seat Code	Seat Width (WW)	Seat Depth (DD)	Seat Configuration (STC) and Access Opening				
<input type="checkbox"/> Z36022M		Molded	22"	<input type="checkbox"/> ZMS 23 18 FRO	23"	18"	<input type="checkbox"/> Front				
<input type="checkbox"/> Z36024M		Molded	24"	<input type="checkbox"/> ZMS 23 18 FRO	23"	18"	<input type="checkbox"/> Front				
<input type="checkbox"/> Z36022C		Custom	22"	<input type="checkbox"/> ZCS ____	____"	____"	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36024C		Custom	24"	<input type="checkbox"/> ZCS ____	____"	____"	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36026C		Custom	26"	<input type="checkbox"/> ZCS ____	____"	____"	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36028C		Custom	28"	<input type="checkbox"/> ZCS ____	____"	____"	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear
<input type="checkbox"/> Z36030C		Custom	30"	<input type="checkbox"/> ZCS ____	____"	____"	<input type="checkbox"/> Front	<input type="checkbox"/> Bridge	<input type="checkbox"/> Left	<input type="checkbox"/> Right	<input type="checkbox"/> Rear

NOTES:

- 1) For models with custom seat types, please specify seat width, depth and seat configuration to designate if the access opening is to be in the Front, Left, Right or Rear. Select Bridge if an access opening is not required. Use this information to create the seat code for the required seat.
- 2) Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for all custom seats.
- 3) Custom seat widths can be 1" or 2" wider than the frame width to provide additional thigh support without increasing overall width of the chair.
- 4) Custom seat depths must be between 16"D and 20"D.
- 5) The Ischial Pelvic Alignment System (IPAS) will allow 1" fore/aft adjustment of the aperture location and seat depth on custom seats. On the molded 23"W x 18"D seat, the IPAS will allow the aperture location to be adjusted up to 2" forward with a corresponding seat depth range from 18"D to 20"D.

USER WEIGHT RANGE

- ☐ Z36WR450 WASS User Weight Range
275 – 450 lbs.

- ☐ Z36WR600 WASS User Weight Range
400 – 600 lbs

BACKREST SELECTION

- ☐ Z563T Multi-Position Back with adjustable-tension vinyl upholstery 20" back height, with 1" fore/aft center of gravity adjustment, two back angles (90°/97°) and flip-up armrests. (Z710HD Not Included)
- ☐ Z710HD Locking Armrest Kit - HD
- ☐ Z5S2 Symphony Backrest - Recessed Solid Back
FRAME WIDTH: ☐ 20" ☐ 22"
- ☐ ZCUSTBK Custom Back Frame Construction

NOTES

SEE OTHER SIDE FOR FOOT & LEG SUPPORT OPTIONS, AND OTHER OPTIONS & ACCESSORIES ►

FOOT AND LEG SUPPORT OPTIONS Z709N ADJUSTABLE FOOTREST INCLUDED WITH Z360

<input type="checkbox"/> Z709N Adjustable Footrest HD 16" – 20.5" (PAIR)	<input type="checkbox"/> Z7ELRL Elevating Legrest – Left INCLUDES RECESSED CALF PAD AND Z7ADMFP
<input type="checkbox"/> Z7FREXT Footrest Extension 18" – 22.5"(PAIR)	<input type="checkbox"/> Z7ELRR Elevating Legrest – Right INCLUDES RECESSED CALF PAD AND Z7ADMFP
<input type="checkbox"/> Z7ADMFP Adjustable Angle/Depth Footplates (PAIR)	<input type="checkbox"/> Z7ELROFL Offset Elevating Legrest – Left (FOR 19" D SEATS) INCLUDES RECESSED CALF PAD AND Z7ADMFP
<input type="checkbox"/> Z712 H-Strap Calf Support HD – 22"W	<input type="checkbox"/> Z7ELROFR Offset Elevating Legrest – Right (FOR 19" D SEATS) INCLUDES RECESSED CALF PAD AND Z7ADMFP
<input type="checkbox"/> Z708ADJ Adjustable Heel Loops (PAIR)	<input type="checkbox"/> ZCUSTFT Custom Foot/Leg Support Construction NOTES
<input type="checkbox"/> Z7OFL Offset Footrest Receiver – Left (MOVES FOOTREST Laterally 2")	
<input type="checkbox"/> Z7OFR Offset Footrest Receiver – Right (MOVES FOOTREST Laterally 2")	

OTHER OPTIONS AND ACCESSORIES

<input type="checkbox"/> Z678 Directional Caster Lock FOR 6" CASTERS	<input type="checkbox"/> Z740AP2L Lateral Support – Left WITH LARGE DEPTH-ADJUSTABLE PAD
<input type="checkbox"/> Z7HCO Headrest Clamp w/ 2" Side-To-Side Offset Adjustment	<input type="checkbox"/> Z740AP2R Lateral Support – Right WITH LARGE DEPTH-ADJUSTABLE PAD
<input type="checkbox"/> Z7HPL Headrest Pad – Large (12" x 6")	<input type="checkbox"/> Z744L Lateral Mounting Hardware – Left REQUIRED FOR LATERAL SUPPORTS
<input type="checkbox"/> Z713 Chest Positioning Strap HD 1-PIECE WITH VELCRO CLOSURE	<input type="checkbox"/> Z744R Lateral Mounting Hardware – Right REQUIRED FOR LATERAL SUPPORTS
<input type="checkbox"/> Z714 Pelvic Positioning Strap HD 2-PIECES WITH FASTEX BUCKLE	<input type="checkbox"/> Z746AP2L Leg Adductor – Left WITH LARGE* DEPTH-ADJUSTABLE PAD
<input type="checkbox"/> Z715L Molded Flat Armpad – Left (14"L x 3.5"W)	<input type="checkbox"/> Z746AP2R Leg Adductor – Right WITH LARGE* DEPTH-ADJUSTABLE PAD
<input type="checkbox"/> Z715R Molded Flat Armpad – Right (14"L x 3.5"W)	<input type="checkbox"/> Z736CXL Chest Strap HD – Bodypoint Aeromesh
<input type="checkbox"/> Z716L Molded Arm Trough – Left	<input type="checkbox"/> Z736PXL Pelvic Strap HD Bodypoint Aeromesh
<input type="checkbox"/> Z716R Molded Arm Trough – Right	<input type="checkbox"/> Z742 Lateral / Adductor 1.5" Extension Plate ea. QUANTITY _____
<input type="checkbox"/> Z717L Pivoting Arm Mount – Left (for Z715L/Z716L)	
<input type="checkbox"/> Z717R Pivoting Arm Mount – Right (for Z715R/Z716R)	
<input type="checkbox"/> ZSPCP1 Commode Pan Spacer Kit 1" Lower	
<input type="checkbox"/> ZSPCP2 Commode Pan Spacer Kit 2" Lower	
<input type="checkbox"/> ZSPCP3 Commode Pan Spacer Kit 3" Lower	
<input type="checkbox"/> ZSPCP4 Commode Pan Spacer Kit 4" Lower	

* LARGE DEPTH-ADJUSTABLE PAD SIZE: 6.75" X 4.3"

CREDIT & OMISSIONS

<input type="checkbox"/> Z905 Omit Commode Pan and Mounting Brackets	<input type="checkbox"/> Z908 Omit Adjustable Footrests (Z709N)
--	---

COMMENTS