

Thigh supports


The key to making the body support self-adjust for you



- SureHands® offers a variety of thigh supports making it possible for many people to be able to use the SureHands® body support.
- Thigh supports fit securely under the thigh and distribute a person's body weight and pressure uniformly. They also ensure that the body maintains posture and is in balance during lifting.
- Thigh supports are made of a maintenance friendly rust-resistant material.
- The patented thigh support fix™ designed for extra comfort.

Thigh supports

● Versions and dimensions



Art. n° 1074 S-M-L-XL (slipproof)
Art. n° 1077 S-M-XL



Art. n° 1075 M (slipproof)
Art. n° 1078 M



Art. n° 1076 S-M (slipproof)
Art. n° 1079 S-M

Art. n°	Description	Color	A [inch/ mm]	B [inch/ mm]
1074S	Standard slipproof small	black	6.30 / 160	5.12 / 130
1074M	Standard slipproof medium	black	8.66 / 220	7.09 / 180
1074L	Standard slipproof large	black	11.22 / 285	9.84 / 250
1074XL	Standard slipproof extra-large	black	12.2 / 310	9.84 / 250
1075M	Contour slipproof medium	black	8.66 / 220	5.91 / 150
1076S	Muscular dystrophy slipproof small	black	6.30 / 160	5.12 / 130
1076M	Muscular dystrophy slipproof medium	black	8.66 / 220	7.09 / 180
1077S	Standard small	grey	6.3 / 160	5.12 / 130
1077M	Standard medium	grey	8.66 / 220	7.09 / 180
1077XL	Standard extra-large	grey	12.2 / 310	9.84 / 250
1078M	Contour medium	grey	8.66 / 220	5.91 / 150
1079S	Muscular dystrophy standard small	grey	6.3 / 160	5.12 / 130
1079M	Muscular dystrophy standard medium	grey	8.66 / 220	7.09 / 180

