

# Pacer dimensions

| Key user dimensions (cm)                                                                                                                                                                                                                                                          |                                  | K610 mini | K620 small | K630 medium | K640 large | K650 XL |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------|-----------|------------|-------------|------------|---------|
| Elbow height                                                                                                                                                                                                                                                                      |                                  | 39–52     | 48–74      | 61–89       | 81–119     | 86–124  |
|  <b>Key user dimension: elbow height</b><br>Measure the vertical distance from the bent <b>elbow to the floor</b> while the user is standing upright.<br>Choose the size that allows for growth. |                                  |           |            |             |            |         |
| Item dimensions (cm)                                                                                                                                                                                                                                                              |                                  | K610 mini | K620 small | K630 medium | K640 large | K650 XL |
| Floor to top of arm prompt                                                                                                                                                                                                                                                        |                                  | 39–52     | 48–74      | 61–89       | 81–119     | 86–124  |
| Floor to top of MPS (not tilted)*                                                                                                                                                                                                                                                 |                                  | N/A       | 25–51      | 24–60       | 43–89      | 50–95   |
| Overall width                                                                                                                                                                                                                                                                     | Standard base:                   | 52        | 58         | 66          | 71         | 80      |
|                                                                                                                                                                                                                                                                                   | Utility base:                    | N/A       | N/A        | 76          | 81         | N/A     |
|                                                                                                                                                                                                                                                                                   | Narrow treadmill/stability base: | N/A       | N/A        | 89          | 89         | N/A     |
|                                                                                                                                                                                                                                                                                   | Wide treadmill/stability base:   | N/A       | N/A        | 102         | 102        | N/A     |
| Overall length                                                                                                                                                                                                                                                                    | Standard base:                   | 57        | 66         | 76          | 91         | 101     |
|                                                                                                                                                                                                                                                                                   | Utility base:                    | N/A       | N/A        | 91          | 103        | N/A     |
|                                                                                                                                                                                                                                                                                   | Treadmill/stability bases:       | N/A       | N/A        | 104         | 104        | N/A     |
| Frame height                                                                                                                                                                                                                                                                      | Standard & utility bases:        | 32        | 41–53      | 55–69       | 72–98      | 79–104  |
|                                                                                                                                                                                                                                                                                   | Treadmill/stability bases:       | N/A       | 58–89      | 60–93       | 76–119     | N/A     |
| Base height (without upper frame)                                                                                                                                                                                                                                                 | Standard & utility bases:        | N/A       | 28         | 37          | 41         | 47      |
|                                                                                                                                                                                                                                                                                   | Treadmill/stability bases:       | N/A       | N/A        | 53          | 53         | N/A     |
| Frame weight (kg)                                                                                                                                                                                                                                                                 | Standard base:                   | 4.3       | 4.7        | 6.8         | 7          | 10      |
|                                                                                                                                                                                                                                                                                   | Utility base:                    | N/A       | N/A        | 8.4         | 8.7        | N/A     |
|                                                                                                                                                                                                                                                                                   | Narrow treadmill/stability base: | N/A       | N/A        | 12.5        | 12.5       | N/A     |
|                                                                                                                                                                                                                                                                                   | Wide treadmill/stability base:   | N/A       | N/A        | 12.7        | 12.7       | N/A     |
|                                                                                                                                                                                                                                                                                   | Standard upper:                  | N/A       | 2.7        | 3.4         | 4.3        | 5       |
|                                                                                                                                                                                                                                                                                   | Dynamic upper:                   | N/A       | 3          | 5.2         | 6          | 6.8     |
| Dynamic upper movement                                                                                                                                                                                                                                                            | Vertical:                        | N/A       | 3.5        | 5           | 5          | 5       |
|                                                                                                                                                                                                                                                                                   | Horizontal:                      | N/A       | 4          | 5           | 5          | 5       |
| Maximum treadmill width                                                                                                                                                                                                                                                           | Narrow treadmill/stability base: | N/A       | N/A        | 74          | 74         | N/A     |
|                                                                                                                                                                                                                                                                                   | Wide treadmill/stability base:   | N/A       | N/A        | 86          | 86         | N/A     |
| Maximum treadmill height                                                                                                                                                                                                                                                          | Treadmill/stability bases:       | N/A       | N/A        | 28          | 28         | N/A     |
| Maximum working load (kg)                                                                                                                                                                                                                                                         |                                  | 22.6      | 34         | 68          | 91         | 113.4   |

\*This measurement is taken at the lowest frame and MPS height, and the highest frame and MPS height.

## Pelvic support

User dimensions (cm)


Indicated by the red arrow below

Key user dimension:

half of thigh circumference


measured at the groin.

| K556 small | K557 med. | K558 large |
|------------|-----------|------------|
| 10–20      | 15–25     | 20–36      |



| Product dimensions (cm) |            |           |            |
|-------------------------|------------|-----------|------------|
|                         | K556 small | K557 med. | K558 large |
| <b>a</b>                | 17         | 22        | 30         |
| <b>b</b>                | 30         | 41        | 51         |
| <b>c</b>                | 18         | 23        | 30         |

## Hip positioner

| Dimensions (cm) – without pad                                                       |            |     |
|-------------------------------------------------------------------------------------|------------|-----|
|  |            |     |
| K513 small                                                                          | K533 large |     |
| <b>a</b> Crotch length                                                              | 6          | 9   |
| <b>b</b> Width                                                                      | 3.5        | 4.5 |
| <b>c</b> Length                                                                     | 23         | 28  |

# Retail Order Form

April 2019

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.

\* Required



Mini and small sizes available Spring 2019.

|                                                   |                                          | NEW K610 mini                 | NEW K620 small                | K630 medium                   | K640 large | K650 XL                      |      |      |
|---------------------------------------------------|------------------------------------------|-------------------------------|-------------------------------|-------------------------------|------------|------------------------------|------|------|
| * Color                                           | CHOOSE                                   | Red <input type="checkbox"/>  | Red <input type="checkbox"/>  | Red                           | Red        | Red <input type="checkbox"/> |      |      |
|                                                   |                                          | Blue <input type="checkbox"/> | Blue                          | Blue                          | Blue       | Blue                         |      |      |
|                                                   |                                          | Gray <input type="checkbox"/> | Gray                          | Gray                          | Gray       | Gray                         |      |      |
|                                                   |                                          | Lime <input type="checkbox"/> | Lime                          | Lime                          | Lime       | Lime                         |      |      |
|                                                   |                                          | Pink <input type="checkbox"/> | Pink                          | Pink                          | Pink       | Pink                         |      |      |
| * Upper frame                                     | CHOOSE                                   | Dynamic                       | N/A                           | K621*                         | K631       | K641                         | K651 |      |
| * K621 is available as of May 2019                |                                          | Standard                      | N/A                           | K622                          | K632       | K642                         | K652 |      |
| * Base                                            | Standard base                            | CHOOSE                        | With odometer                 | K613 <input type="checkbox"/> | K623       | K633                         | K643 | K653 |
|                                                   |                                          | Without odometer              | K614 <input type="checkbox"/> | K624                          | K634       | K644                         | K654 |      |
|                                                   | Utility base (larger wheels)             | CHOOSE                        | With odometer                 | N/A                           | N/A        | K635                         | K645 | N/A  |
|                                                   |                                          | Without odometer              | N/A                           | N/A                           | K636       | K646                         | N/A  |      |
|                                                   | Treadmill/stability base (comes in gray) | CHOOSE                        | Standard                      | N/A                           | K676       | K676                         | K676 | N/A  |
|                                                   |                                          | Wide                          | N/A                           | K678                          | K678       | K678                         | K678 | N/A  |
| Arm supports (pair)                               | CHOOSE                                   | Small arm prompts             | K512                          | K512                          | K512       | K512                         | K512 |      |
|                                                   |                                          | Large arm prompts             | N/A                           | N/A                           | K532       | K532                         | K532 |      |
|                                                   |                                          | Arm platforms                 | K697                          | K697                          | K697       | K697                         | K697 |      |
|                                                   |                                          | Arm platforms w/ handgrips    | K696                          | K696                          | K696       | K696                         | K696 |      |
| Handbrake for utility base (requires arm prompts) |                                          | N/A                           | N/A                           | K687                          | K688       | N/A                          |      |      |
| Hand loops (pair)                                 |                                          | K579                          | K579                          | K579                          | K579       | K579                         |      |      |
| Chest prompt                                      | CHOOSE                                   | Small                         | K514                          | K514                          | K514       | N/A                          | N/A  |      |
|                                                   |                                          | Medium                        | K524                          | K524                          | K524       | K524                         | K524 |      |
|                                                   |                                          | Large                         | N/A                           | N/A                           | K534       | K534                         | K534 |      |
| Multi-position saddle                             | CHOOSE                                   | New! Small                    | N/A                           | K672                          | N/A        | N/A                          | N/A  |      |
|                                                   |                                          | Large                         | N/A                           | N/A                           | K674       | K674                         | K674 |      |
| Pelvic support (includes handholds)               | CHOOSE                                   | Small                         | K556                          | K556                          | K556       | N/A                          | N/A  |      |
|                                                   |                                          | Medium                        | N/A                           | K557                          | K557       | K557                         | N/A  |      |
|                                                   |                                          | Large                         | N/A                           | N/A                           | K558       | K558                         | K558 |      |
| Hip positioner (includes handholds)               | CHOOSE                                   | Small                         | K513                          | K513                          | K513       | K513                         | K513 |      |
|                                                   |                                          | Large                         | N/A                           | N/A                           | K533       | K533                         | K533 |      |
| Hip positioner pad                                | CHOOSE                                   | Small                         | K541                          | K541                          | K541       | K541                         | K541 |      |
|                                                   |                                          | Large                         | N/A                           | N/A                           | K542       | K542                         | K542 |      |
| Thigh prompts (pair)                              | CHOOSE                                   | Small                         | N/A                           | K515                          | K515       | K515                         | N/A  |      |
|                                                   |                                          | Large                         | N/A                           | N/A                           | K535       | K535                         | K535 |      |
| Ankle prompts for standard base (pair)            |                                          | K681                          | K682                          | K683                          | K684       | K684                         |      |      |
| Ankle prompts for utility base (pair)             |                                          | N/A                           | N/A                           | K693                          | K694       | N/A                          |      |      |
| Ankle prompts for treadmill/stability base (pair) |                                          | K664                          | K664                          | K664                          | K664       | N/A                          |      |      |
| Additional handholds (pair)                       |                                          | K519                          | K519                          | K519                          | K519       | K519                         |      |      |
| Communication tray                                |                                          | K518                          | K698                          | K698                          | K698       | K698                         |      |      |
| Attendant guide bar                               |                                          | K589                          | K679                          | K689                          | K689       | N/A                          |      |      |
| Rifton accessories tote                           |                                          | K522                          | K522                          | K522                          | K522       | K522                         |      |      |
| Total                                             |                                          | K610                          | K620                          | K630                          | K640       | K650                         |      |      |