

Cruiser X4 Wheelchair

HCPCS: K0004-1,K0004



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Weighs 34 pounds or less when properly equipped in the basic configuration (Excluding front rigging)
- Easy to propel and transport
- Frame style eliminates seat guides and allows for custom backs and accessories
- Flip-back, removable arm is height adjustable from top of seat rail 10" to 14" and features a front mounted easy-to-use lever release mechanism (Figure A - C)
- Padded armrests provide added patient comfort
- Height adjustable segmented back
- Aluminum and Carbon steel blended frame with combination powder coated and silver vein finish is attractive, chip resistant and easy-to-maintain
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Patented adjustable depth seat 16" and 18" on chair depths (Figure D - F)
- 3 position quick release rear axles mounted on Composite Mag-style wheels with composite hand rims are lightweight and maintenance free (Figure G)
- Polyurethane rear tires and composite casters offer superior performance, smoother ride and are lightweight
- Dual over center cross links provide added rigidity (Figure H)
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Composite footplates with heel loops are durable and lightweight
- Precision sealed wheel bearings through-out ensure long lasting and reliability
- 8" front casters have 3 height adjustments and angle adjustments (Figure I)
- Triple rear axle positions provides easy transition of seat height from standard to hemi level
- Comes with push-to-lock wheel locks

PRODUCTS

Description

CX416ADDA-ELR	CruiserX4, 16" Adj DeskArm, ElvLg, 1/ea
CX416ADDA-SF	CruiserX4, 16" Adj DeskArm, SwgFt, 1/ea
CX416ADFA-ELR	CruiserX4, 16" Adj FullArm, ElvLg, 1/ea
CX416ADFA-SF	CruiserX4, 16" Adj FullArm, SwgFt, 1/ea
CX418ADDA-ELR	CruiserX4, 18" Adj DeskArm, ElvLg, 1/ea
CX418ADDA-SF	CruiserX4, 18" Adj DeskArm, SwgFt, 1/ea
CX418ADFA-ELR	CruiserX4, 18" Adj FullArm, ElvLg, 1/ea
CX418ADFA-SF	CruiserX4, 18" Adj FullArm, SwgFt, 1/ea
CX420ADDA-ELR	CruiserX4, 20" Adj DeskArm, ElvLg, 1/ea
CX420ADDA-SF	CruiserX4, 20" Adj DeskArm, SwgFt, 1/ea
CX420ADFA-ELR	CruiserX4, 20" Adj FullArm, ElvLg, 1/ea
CX420ADFA-SF	CruiserX4, 20" Adj FullArm, SwgFt, 1/ea

SPECIFICATIONS

Description

Back Height	16" - 17.5"
Carton Shipping Weight	46 lbs. 48 lbs. 50 lbs.
Height	36"
Length with Riggings	49"

Customer Service | 866.752.1363

Product Weight without Riggings (16", 18", 20") 32 lbs. | 34 lbs. | 36 lbs.