

Rotating Suction-Cup Grab Bar

HCPCS: E0241



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Installs and removes without tools or professional installation
- Rotates to any angle (Figure A)
- Will not damage property
- Large suction cups provide an extremely strong hold
- Release levers make installing and removing the suction-cup grab bars quick and easy
- Red/green color indicator shows "safe secure" suction or "re-affix" suction (Figure B - C)
- Avoid grout lines for proper placement
- Retail packaged

PRODUCTS

Description

RTL13084	Grab Bar SuctionCup,ADJ,19 3/4",1EA,3/CS
----------	--

SPECIFICATIONS

Description

Carton Shipping Weight	3.5 lbs.
Length	19.75"
Warranty	Limited lifetime
Weight	1.0 lbs.
Grab Bars	Length 19.75"