

# SAFETYSURE® CURVED TRANSFER BOARD

## *Curved Design for Easier Transfers*



Available in either Plastic or Wood

## SafetySure® Curved Transfer Board

The arch shape of the SafetySure® Curved Transfer Board enables a variety of positioning options but also allows for the natural sweeping motion of an assisted transfer. The Curved Transfer Board is ideal for bed, car, wheelchair and toileting transfers. Available in either plastic or wood. The plastic board is made from durable, semi-rigid polyethylene that can be easily wiped down for infection control. The wood board is made from plywood and finished with a lacquered coating. The surface of each board provides a smooth sliding surface and the ends are flared to provide a stable transfer platform. The Curved Transfer Board can be used for assisted and/or unassisted transfers. Weight capacity 250 lbs.



### Product Specifications:

Item #	Material	Width (Center)	Width (End)	Length	Product Weight	Weight Capacity
7005837	Plastic	10.5"	12"	28"	2.6 lbs.	250 lbs.
7005838	Wood	10.5"	12"	28"	3.5 lbs.	250 lbs.