

Room to grow

Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 15"
(with footboard lift 11")



Large chair at its largest dimensions

- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 20"



Large standard base
R860 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
	Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.
	Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.
Item dimensions (inches)	R860 large Standard base
Frame width	short legs: 26¾
	long legs: 26¾
	short legs w/ casters: 27¾
	long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½
	long legs: 20–24
	short legs w/ casters: 16½–20½
	long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	15–20
	11–16
Seat width with hip guides without hip guides	11–14
	17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lbs)	250



Large hi/lo base
R870 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
	Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.
	Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.
Item dimensions (inches)	R870 large Hi/lo base
Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	15–20
	11–16
Seat width with hip guides without hip guides	11–14
	17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lbs)	250



Retail Order Form
July 2018

All Rifton products are made to order and are non-returnable and non-exchangeable.



		Standard R860 large		Hi/lo R870 large			
Required components	Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE Adjustable	R861	<input type="checkbox"/>	R871	<input type="checkbox"/>	
		Adjustable w/ back spring	R862	<input type="checkbox"/>	R872	<input type="checkbox"/>	
	Base	CHOOSE Without spring	R873	<input type="checkbox"/>	R883	<input type="checkbox"/>	
		With spring	R874	<input type="checkbox"/>	N/A		
	Arm supports (pair)	CHOOSE Armrests	R868	<input type="checkbox"/>	R868	<input type="checkbox"/>	
		Forearm prompts	R879	<input type="checkbox"/>	R879	<input type="checkbox"/>	
	Adjustable legs	CHOOSE	Short	R808	<input type="checkbox"/>	N/A	
			Long	R894	<input type="checkbox"/>	N/A	
			Short with casters	R895	<input type="checkbox"/>	N/A	
			Long with casters	R896	<input type="checkbox"/>	N/A	
	Pads	CHOOSE	Tan	R875	<input type="checkbox"/>	R875	<input type="checkbox"/>
			Blue	R876	<input type="checkbox"/>	R876	<input type="checkbox"/>
			Red	R877	<input type="checkbox"/>	R877	<input type="checkbox"/>
			Pink	R977	<input type="checkbox"/>	R977	<input type="checkbox"/>
Purple			R978	<input type="checkbox"/>	R978	<input type="checkbox"/>	
Green			R979	<input type="checkbox"/>	R979	<input type="checkbox"/>	
Optional accessories	Headrest	CHOOSE	Contoured	R897	<input type="checkbox"/>	R897	<input type="checkbox"/>
			Flat	R898	<input type="checkbox"/>	R898	<input type="checkbox"/>
			Adjustable winged	R881	<input type="checkbox"/>	R881	<input type="checkbox"/>
	Laterals (pair)	CHOOSE	Small	R828	<input type="checkbox"/>	R828	<input type="checkbox"/>
			Large	R838	<input type="checkbox"/>	R838	<input type="checkbox"/>
			Large with strap	R839	<input type="checkbox"/>	R839	<input type="checkbox"/>
	Push handles (recommended for all chairs with wheels)		R878	<input type="checkbox"/>	included		
	Tray (requires armrests)		R669	<input type="checkbox"/>	R669	<input type="checkbox"/>	
	Handhold		R649	<input type="checkbox"/>	R649	<input type="checkbox"/>	
	Chest strap (wide)		R807	<input type="checkbox"/>	R807	<input type="checkbox"/>	
	Butterfly harness	CHOOSE	Regular	R889	<input type="checkbox"/>	R889	<input type="checkbox"/>
			Slim cut	R826	<input type="checkbox"/>	R826	<input type="checkbox"/>
	Thigh belt		R813	<input type="checkbox"/>	R813	<input type="checkbox"/>	
	Pelvic harness		R846	<input type="checkbox"/>	R846	<input type="checkbox"/>	
	Hip guides (pair)		R867	<input type="checkbox"/>	R867	<input type="checkbox"/>	
	Abductor	CHOOSE	Small	R884	<input type="checkbox"/>	R884	<input type="checkbox"/>
			Large	R891	<input type="checkbox"/>	R891	<input type="checkbox"/>
	Adductors (pair)		R892	<input type="checkbox"/>	R892	<input type="checkbox"/>	
	Leg prompt		R859	<input type="checkbox"/>	R859	<input type="checkbox"/>	
	Footboard (recommended for chairs with long legs)		R857	<input type="checkbox"/>	included		
	Footboard lift		R818	<input type="checkbox"/>	R818	<input type="checkbox"/>	
	Ankle straps (pair)		R855	<input type="checkbox"/>	R855	<input type="checkbox"/>	
	Sandals (pair)		R964	<input type="checkbox"/>	R964	<input type="checkbox"/>	
	Wedges (pair, requires sandals)		R699	<input type="checkbox"/>	R699	<input type="checkbox"/>	
	Lumbar and seat support kit		R888	<input type="checkbox"/>	R888	<input type="checkbox"/>	
	Backrest filler pad		R856	<input type="checkbox"/>	R856	<input type="checkbox"/>	
	Whitmyer headrest adapter plate		R804	<input type="checkbox"/>	R804	<input type="checkbox"/>	
	Total		R860		R870		