


How to buy

- 1. Measure** your client to determine the size you need.
(See dimension chart for help.)
- 2. Choose** the accessories you need.



Retail Order Form December 2018

All Rifton products are made to order and are non-returnable and non-exchangeable.

User dimensions (cm)	R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
Height	91 – 119	109 – 137	124 – 152	140 – 173	157+
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth, and height are adequate for each user, and that the user's weight does not exceed the maximum weight recommended.					
Item dimensions (cm)	R310	R320	R330	R340	R350
Seat height	23,25,28	27,29,32	32,34,37	37,39,42	42,44,47
Seat depth	24	28	32.5	37.5	42
Seat width	27	31	36	41.5	46.5
Backrest height	24	28	32	36.5	41
Armrest height above the seat	14	15	18	21	22.5
Armrest length	17	19	22	24.5	28
Max. working load (kg)	45	55	80	100	113



	R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
* Compass chair	R311 <input type="checkbox"/>	R321 <input type="checkbox"/>	R331 <input type="checkbox"/>	R341 <input type="checkbox"/>	R351 <input type="checkbox"/>
Seat belt	R312 <input type="checkbox"/>	R312 <input type="checkbox"/>	R332 <input type="checkbox"/>	R342 <input type="checkbox"/>	R342 <input type="checkbox"/>
Stability feet	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>
Total	R310	R320	R330	R340	R350

* Required