

Drive Full Electric Low Height Bed

HCPCS: E0297



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F

- Low Deck height with flexibility to configure "roll at any height" or "roll at low height only" castor location
- -Roll at any Height: 14.81" (376mm) to 22.83" (580 mm) deck height
- -Roll at Low Height: 11.81 (300 mm) to 19.83" (504mm) deck height
- Quiet, convenient, full electric operation provides smooth position change w/easy-to-use Hand Pendant control of Head/Foot Section positioning and overall height of the bed frame for maximum comfort and versatility for the patient or caregiver
- Height-adjustment and Head/Foot deck motors come pre-installed on the foot section no assembly required
- In the event of power failure, top mounted Hi/Lo crank accessibility and manual crank motor override of sleep surface adjustments
- Easy-to-transport and set up, Bed ships in two cartons
- Channel frame construction provides superior strength and reduced weight
- Reinforced frame resists twisting and bending
- Long lasting Zinc-coated spring deck for Patient comfort
- Heavy-duty bed ends with durable HD Polyethylene Molded Bed End covers are attractive, easy-to-maintain
- Headboard is taller than footboard to blend with home décor
- Labels and color-coded springs help ensure accurate installation of side rails for patient safety
- High performance workhorse muscle of the 3-motor power plant is an industry standard
- Half rails and full rails available in packages (PKG) or sold separately
- Motor is UL Approved, Bed is manufactured in ISO9001 & ISO13485 facility
- Satisfies applicable FDA entrapment guidelines
- Extension kit available (model #15005EXTKIT-L) extends length from 80" to 84"
- Patient Weight Capacity 350 lbs. (158kg). Total Weight Capacity 450 lbs. (204kg). Includes occupant, mattress, rails, trapeze, accessories, personal affects, etc.

PRODUCTS

Description

15005LP	Full Elec LowBed 1/cv
15005LP-FR	Full Elec LowBed,Full Rail Pkg
15005LP-FRPKG	Full Elec LowBed,FRail IMatt Pkg
15005LP-FRPKG-FM	Full Elec LowBed,FRail FMatt Pkg
15005LP-FRPKG-TS	Full Elec LowBed,FRail TSMatt Pkg
15005LP-HR	Full Elec LowBed Half Rail Pkg
15005LP-HRPKG	Full Elec LowBed HRail IMatt Pkg
15005LP-HRPKG-TS	Full Elec LowBed HRail TSMatt Pkg

SPECIFICATIONS

Description

Carton Shipping Weight	184 lbs
Roll at Low Height	11.81 (300mm) to 19.83 (504mm) deck height

SPECIFICATIONS**Description**

Roll at any Height	14.81" (376mm) to 22.83 (580mm) deck height
Unit Dimensions	36" (.914m) (W) x 88" (2.24m) (L)
Warranty	Original Purchaser. Limited Lifetime on welds, 5 year on frame, motor, pendant, 1 year on all other parts and components. During the warranty period, defective items will be repaired or replaced at Drive's option at no charge
Weight	169 lbs