

INVACARE®

# AQUATEC® COMMODES

## OCEAN ERGO VIP TILT IN SPACE

---



---

Customer insight-led design to enhance independence, safety and flexibility. The Aquatec® Ocean Ergo Family of shower chair commodes have been designed with user comfort, independence and dignity in mind. Based on the well-established Aquatec® Ocean range, the Ocean Ergo family comes packed with a range of new features designed to make personal care activities easier for both the user and caregiver

The Ocean Ergo VIP boasts a more ergonomic design and features low force center of gravity tilt with a range of -5 to 40 degrees of tilt. The new Ergo VIP provides a smooth transition through the tilt range without the use of a gas cylinder and reduces with overall foot print of the chair while providing maximum stability in all positions. The new frame design also enables Ergo VIP to get closer to the wall to better position the shower chair commode over a toilet for daily hygiene needs.

---

## FEATURES

- Ergonomic seat plate with a 5° posterior angle to support natural toileting position
- Piston free tilting mechanism reduces effort required to position centre of gravity
- Enhanced frame design for greater over toilet positioning and forward stability
- Collection pan with integrated measuring scale included
- Tool-free seat height adjustment range of 5" to fit over most toilets
- Rust proof, stainless steel frame
- Machine washable adjustable tension backrest

---

## ACCESSORIES



1. Soft Seat with Hygiene Recess
2. Universal soft seat - Ocean Ergo Series
3. Front or Rear Anti-tipper
4. Adjustable Leg-rest Assembly (available in Right or Left)
5. Hygiene Backrest Support

---

## SPECIFICATIONS

MODEL	ERGO VIP
Length	39"
Width	22"
Height	40.5" - 45.5"
Product Weight	51.8 lbs
Width Between Arms/Handles	18"
Seat Width	19"
Back Width	19"
Seat Depth	19"
Seat Angle	5°
Seat to Footrest	12.5" - 18"
Height Seat to Floor	18.5" - 23"
Tilt Range	-5 - 40°
Capacity	330 lbs

---

## WARRANTY

Aquatec® OceanErgo Series shower commodes are covered under a two-year limited warranty against manufacturing defects.

**Manufacturer's warranty does not cover normal wear.**

All Aquatec® OceanErgo shower commode accessories are covered for a period of one year.