

CHILD RESTRAINT SYSTEM

CARROT 3

BOOSTER SEAT




- **Adaptive Positioning**
High density foam inserts for optimal positioning and maximum comfort
- **Grows With Your Child**
Able to accommodate growth from 79-165 lbs / 55"-69"
- **"Free Angle" Recline up to 180°**
Conforms to the vehicle's seat
- **LATCH Equipped**
Lower anchors and tethers for children for easy and safe installation
- **Offers a Large Configuration**
To meet the needs of children and adults who require back height support and a wide seat

Shoulder Protector Wings



Gives lateral shoulder support and protection
(Width Adjustable)

Wide Seat



The Carrot 3 Booster Seat is for users who require a wide seat (16.5")

4-Point Positioning H-Harness



To maintain midline positioning of the trunk for increased postural control


Booster Kit



Option to transform an existing Carrot 3 into a booster in 3 easy steps


Lightweight • Compact Folding • Convenient

Carrot 3 Booster Kit
(Option to transform Carrot 3 Seat to Booster Seat)



- 1 **Booster Seat Base**
- 2 **4-Point Positioning H-Harness**
- 3 **H-Harness Attaching Bracket**
- 4 **Backrest Attaching Metal Rods**

Carrot 3 Booster Kit
Required Extensions



- 1 **Shoulder Protector Wings**
- 2 **2" Backrest Extension**


Accessories



Upholstery Colors



Carrot 3 Specifications with Extensions



*Height is measured from the top of the headrest to the seat.

Booster Seat Standard Specifications

	Seat Width	Seat Depth	Backrest Height	Head Support Height	Recline Adjustment	Product Weight	Product Overall Width
Measurements:							
INCHES	16.5	17	28 - 33	8.5	90° - 180°	21 lbs	20 - 22
	User Weight Range	User Height Range	Shoulder Protector Wings Height	Backrest Height	Backrest Extensions Height	Shoulder Wing Width Adjustment	
INCHES	79 lbs - 165 lbs	55 - 69	6	11.8	2 or 4	15.25 - 19.25	

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately. Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.