

## DV8 Steerable Aluminum Knee Walker

HCPCS: E0118



FIGURE A



FIGURE B




FIGURE C



FIGURE D



FIGURE E




FIGURE F



FIGURE G

- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- Deluxe dual-braking system
- Sturdy dual frame adds stability and durability
- New basket allows the user to transport personal items (Figure D)
- Contoured leg pad provides extra balance
- New leg pad cover provides added comfort (Figure E)
- New tiller lock provides added safety and security
- Tiller can be easily folded for storage and transportation
- Foldable steering column allows for easy storage
- Tool-free height adjustment of seat and handle (Figure A - B)
- 8" casters are ideal for indoor/outdoor use
- Silver vein finish (Figure C) is durable and easy-to-maintain

### PRODUCTS

### Description

790	Steerable Knee Walker,1/ea
-----	----------------------------

### SPECIFICATIONS

### Description

Carton Shipping Weight	26 lbs.
Handle Height	33" - 37.5"
Knee Pad Height	17.5" - 21.5"
Unit Dimensions	17" (W) x 32" (D)
Warranty	Limited Lifetime, Lifetime Warranty on Brake Cable
Weight	19 lbs.
Weight Capacity	300 lbs

### PARTS

### Description

9503H79015-N	KneePadHeightKnob790New,1ea
9503H79016	Left Wheel Lock,790,1ea
9503H79017	Right Wheel Lock,790,1ea
9503H79018	Caster Complete Assy,790,1ea
9503H79019	Caster Axle,790,1ea
9503H79020	Axle Set Screws,790,4/pkg
9503H79021	Wheel Bearings,790,PhoenixHD3,1/pr
9503H79001	Handle Bar,790,1ea
9503H79002	Hand Grips,790,1 Pair

## PARTS

## Description

9503H79003	Left Brake Complete,790,1ea
9503H79004	Right Brake Complete,790,1ea
9503H79005	Handle Bar Column,790,1ea
9503H79006	Handle Bar Height Adj Lock,790,1ea
9503H79007	Steering Column,790,1ea
9503H79008	Brake Cable Guide,790,1ea
9503H79002B	
9503H79011	Steering Fork,790,1ea
9503H79012	Column To Fork Locking Bolt,790,1ea
9503H79013	Fork Bearings,790,1pair
9503H79023	Brake Cable & Hardware Only for 790, 1ea