



JAY GO Cushion

Effective July 2018




How to Order:

Part #
☐ **JGO1818** **JAY GO Cushions**
JAY GO Cushion 18" x 18"

Each JAY GO part number is simply based upon the size you are ordering.
You may also reference the size matrix on page 1 of this price list.

Part # Example



JAY GO Cushion Cover

The JAY GO cover is designed using dartex™ fabric which is known to be both stretchable and water resistant. This cover's sizing and sewing pattern is designed to reduce shear and peak pressure buildup. The material itself provides stretch which is an uncommon property in most incontinent products.



Figure A

Please visit www.dartexcoatings.com for additional information on the use of dartex coatings



JAY GO Cushion Foam

The JAY GO foam base is designed to contain an accommodating seating surface in a thicker, softer dual layered foam base. The optimal combination of two different foam densities provides skin protection.



Cross section view

JAY GO Size Matrix												
	WIDTH											
		14	15	16	17	18	19	20	21	22	23	24
Depth	14											
	15											
	16											
	17											
	18											
	19											
	20											
	21											
	22											
	23											
	24											

NO Modified sizes are available

☐ Size is available

☐ Size is not available

JAY GO Cushion Sealed Zipper

The JAY GO zipper and cover provide protection to the foam base against odor because the location is low and away from fluid. The water resistant material, combined with the Aquaguard™ zipper and anti wicking thread provides a thorough protection system for a cushion base.



Figure B

Optimized Aquaguard Zipper location

The JAY GO Cushion is an excellent complement to the JAY GO Back.

Please see the **JAY GO Back** order form (part # 119284) for additional information.



Each JAY GO cushion includes an outer cover, dual layered foam base, plastic bag with label and Owner's Manual


JAY GO Cushion

Effective July 2018

Part #	Description
☐ JGO1414	JAY GO Cushion 14" x 14"
☐ JGO1415	JAY GO Cushion 14" x 15"
☐ JGO1416	JAY GO Cushion 14" x 16"
☐ JGO1417	JAY GO Cushion 14" x 17"
☐ JGO1418	JAY GO Cushion 14" x 18"
☐ JGO1419	JAY GO Cushion 14" x 19"
☐ JGO1420	JAY GO Cushion 14" x 20"
☐ JGO1514	JAY GO Cushion 15" x 14"
☐ JGO1515	JAY GO Cushion 15" x 15"
☐ JGO1516	JAY GO Cushion 15" x 16"
☐ JGO1517	JAY GO Cushion 15" x 17"
☐ JGO1518	JAY GO Cushion 15" x 18"
☐ JGO1519	JAY GO Cushion 15" x 19"
☐ JGO1520	JAY GO Cushion 15" x 20"
☐ JGO1521	JAY GO Cushion 15" x 21"
☐ JGO1614	JAY GO Cushion 16" x 14"
☐ JGO1615	JAY GO Cushion 16" x 15"
☐ JGO1616	JAY GO Cushion 16" x 16"
☐ JGO1617	JAY GO Cushion 16" x 17"
☐ JGO1618	JAY GO Cushion 16" x 18"
☐ JGO1619	JAY GO Cushion 16" x 19"
☐ JGO1620	JAY GO Cushion 16" x 20"
☐ JGO1621	JAY GO Cushion 16" x 21"
☐ JGO1622	JAY GO Cushion 16" x 22"
☐ JGO1714	JAY GO Cushion 17" x 14"
☐ JGO1715	JAY GO Cushion 17" x 15"
☐ JGO1716	JAY GO Cushion 17" x 16"
☐ JGO1717	JAY GO Cushion 17" x 17"
☐ JGO1718	JAY GO Cushion 17" x 18"
☐ JGO1719	JAY GO Cushion 17" x 19"
☐ JGO1720	JAY GO Cushion 17" x 20"
☐ JGO1721	JAY GO Cushion 17" x 21"
☐ JGO1722	JAY GO Cushion 17" x 22"
☐ JGO1723	JAY GO Cushion 17" x 23"
☐ JGO1814	JAY GO Cushion 18" x 14"
☐ JGO1815	JAY GO Cushion 18" x 15"
☐ JGO1816	JAY GO Cushion 18" x 16"
☐ JGO1817	JAY GO Cushion 18" x 17"
☐ JGO1818	JAY GO Cushion 18" x 18"
☐ JGO1819	JAY GO Cushion 18" x 19"
☐ JGO1820	JAY GO Cushion 18" x 20"
☐ JGO1821	JAY GO Cushion 18" x 21"
☐ JGO1822	JAY GO Cushion 18" x 22"
☐ JGO1823	JAY GO Cushion 18" x 23"
☐ JGO1824	JAY GO Cushion 18" x 24"
☐ JGO1915	JAY GO Cushion 19" x 15"
☐ JGO1916	JAY GO Cushion 19" x 16"
☐ JGO1917	JAY GO Cushion 19" x 17"
☐ JGO1918	JAY GO Cushion 19" x 18"
☐ JGO1919	JAY GO Cushion 19" x 19"
☐ JGO1920	JAY GO Cushion 19" x 20"
☐ JGO1921	JAY GO Cushion 19" x 21"
☐ JGO1922	JAY GO Cushion 19" x 22"
☐ JGO1923	JAY GO Cushion 19" x 23"
☐ JGO1924	JAY GO Cushion 19" x 24"
☐ JGO2015	JAY GO Cushion 20" x 15"
☐ JGO2016	JAY GO Cushion 20" x 16"
☐ JGO2017	JAY GO Cushion 20" x 17"
☐ JGO2018	JAY GO Cushion 20" x 18"
☐ JGO2019	JAY GO Cushion 20" x 19"
☐ JGO2020	JAY GO Cushion 20" x 20"
☐ JGO2021	JAY GO Cushion 20" x 21"
☐ JGO2022	JAY GO Cushion 20" x 22"
☐ JGO2023	JAY GO Cushion 20" x 23"
☐ JGO2024	JAY GO Cushion 20" x 24"
☐ JGO2117	JAY GO Cushion 21" x 17"
☐ JGO2118	JAY GO Cushion 21" x 18"
☐ JGO2119	JAY GO Cushion 21" x 19"
☐ JGO2120	JAY GO Cushion 21" x 20"

Part #	Description
☐ JGO2121	JAY GO Cushion 21" x 21"
☐ JGO2122	JAY GO Cushion 21" x 22"
☐ JGO2123	JAY GO Cushion 21" x 23"
☐ JGO2124	JAY GO Cushion 21" x 24"
☐ JGO2217	JAY GO Cushion 22" x 17"
☐ JGO2218	JAY GO Cushion 22" x 18"
☐ JGO2219	JAY GO Cushion 22" x 19"
☐ JGO2220	JAY GO Cushion 22" x 20"
☐ JGO2221	JAY GO Cushion 22" x 21"
☐ JGO2222	JAY GO Cushion 22" x 22"
☐ JGO2223	JAY GO Cushion 22" x 23"
☐ JGO2224	JAY GO Cushion 22" x 24"
☐ JGO2319	JAY GO Cushion 23" x 19"
☐ JGO2320	JAY GO Cushion 23" x 20"
☐ JGO2321	JAY GO Cushion 23" x 21"
☐ JGO2322	JAY GO Cushion 23" x 22"
☐ JGO2323	JAY GO Cushion 23" x 23"
☐ JGO2324	JAY GO Cushion 23" x 24"
☐ JGO2419	JAY GO Cushion 24" x 19"
☐ JGO2420	JAY GO Cushion 24" x 20"
☐ JGO2421	JAY GO Cushion 24" x 21"
☐ JGO2422	JAY GO Cushion 24" x 22"
☐ JGO2423	JAY GO Cushion 24" x 23"
☐ JGO2424	JAY GO Cushion 24" x 24"

JAY GO MODIFICATIONS

MOD	Description	MOD CODE
☐	Custom Embroidery Must specify location and font? Location _____ Font _____ Includes up to 32 characters in one single spot Write specific text in "NOTES" box below	
☐	Notch for Seat Rails 1.25" from bottom of cushion removed 1.25" from side of cushion removed 	NR
☐	Notch for Posts (Back Canes) 1" from side of cushion removed Must specify 1, 2 or 3" from back of cushion  Notch depth _____ No modification to inner or outer cover	NP
☐	Anterior or Posterior Slope Must specify posterior or anterior and height Foam wedge glued underneath cushion, covers sewn to match Anterior or Posterior _____ Height increase _____	WGP or WGA
☐	Leg Length Cut	LLCR/LLCL
☐	Contracture Cut (1" x 45°)	CCT
☐	Leg Trough Cut (1" std)	LTC
☐	Custom Drop Seat	CDS
☐	Medial Support Cut	MSC
☐	Color Cover	CC
☐	Scrotal Cut	SC

Notes:

Please visit SunParts Online for spare parts information

