

2000IPS™

CAR SEAT



ACCESSORIES



Extensor Thrust Wedge



Trunk Support Pads



Buckle Guard



Retainer Clip Guard



Lateral Pads



Bolt-In Hardware Kit

Inspired
by drive



The lightweight and narrow IPS™ Car Seat provides the safety of a traditional car seat while providing a high level of postural support for users with mild positioning needs.

BENEFITS

Adapts to User

FEATURES

- Supports users up to 60" in height and 20–102 lb.
- Three sizes with seat depth extenders (10, 12.5, and 15") grow with user

Positioning Support

- Height-adjustable headwings come standard
- Memory foam back/seat padding reduces sores and fatigue

Built-In Safety

- The IPS Car Seat offers a lighter and narrower alternative to our Spirit Car Seat. Designed to accommodate children with mild positioning needs, the IPS offers the safety of an everyday car seat while providing a higher level of postural support than a traditional car seat.

Versatile & Attractive

- Available in Speedway Gray
- Approved for airline travel, buses, vans, and cars
- Removable and washable upholstery

For more information, call us today!

877-224-0946 | www.inspiredbydrive.com