



FEATURES

- Accommodate users weighing 20 - 102 lb and up to 60" tall
- Durable, machine washable, comfortable padded cover in three color options
- Crash tested to FMVSS #213
- Certified for use in aircraft
- Positioning options available
- Optional seat depth extender allows car seat to grow with the child

INCLUDES

- Height adjustable soft-padded headwings
- Memory foam back/seat padding reduces sores and fatigue
- Standard 5-point harness restraint for better protection than vehicle lap/shoulder belt
- Crotch strap pad
- Lateral positioning pads
- Dual-hook tether strap

UPHOLSTERY



Convertible Pink
(RR)



Roadster Red
(RR)



Speedway Gray
(SG)

ACCESSORIES



Extensor Thrust
Wedge



Trunk Support Pads



Buckle Guard



Retainer Clip Guard

ORDER FORM

FAX TO 905-787-8550 OR SUBMIT ONLINE

Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____

SIZE	Standard Seat Depth	Standard Seat Depth w/Abductor	Medium Seat Depth	Medium Seat Depth w/Abductor	Large Seat Depth	Large Seat Depth w/Abductor
Speedway Gray SG	CSI-2000W-SG	CSI-2000WA-SG	CSI-2000WM-SG	CSI-2000WMA-SG	CSI-2000WL-SG	CSI-2000WLA-SG
Roaster Red RR	CSI-2000W-RR	CSI-2000WA-RR	CSI-2000WM-RR	CSI-2000WMA-RR		
Convertible Pink CP	CSI-2000W-CP	CSI-2000WA-CP				
SELECT TO ITEMIZE		CSI-2000WA-XX-IT	CSI-2000WM-XX-IT	CSI-2000WMA-XX-IT	CSI-2000WL-XX-IT	CSI-2000WLA-XX-IT
Seat Depth Extender			☒ CSI-K2025	☒ CSI-K2025	☒ CSI-K2026	☒ CSI-K2026
Abductor		☒ CSI-K2015		☒ CSI-K2018		☒ CSI-K2019
Weight Capacity	20 - 102 lb	20 - 102 lb	20 - 102 lb	20 - 102 lb	20 - 102 lb	20 - 102 lb
Max. Height Capacity	60"	60"	60"	60"	60"	60"
Harness Height	11.5", 14.5", 17.5", 20.5"					
Back Height	29"	29"	29"	29"	29"	29"
Seat Depth	10"	10"	12.5"	12.5"	15"	15"
Inside Seat Width, Front	15.5"	15.5"	15.5"	15.5"	15.5"	15.5"
Inside Seat Width, Rear	13"	13"	13"	13"	13"	13"
Overall Width	19"	19"	19"	19"	19"	19"
Overall Depth (w/seat rails)	13.5"	16.5"	16"	19"	18.5"	21.5"
Seat-to-Back Angle	88°	88°	88°	88°	88°	88°
Product Weight	16	18	19	22	20	25

ACCESSORIES (For all IPS car seat models)

Block Style Trunk Support Pads
(2 each of 1.5" and 2.5" thick pads)

CSI-A2008-BK

Lateral Pads (set of 4)

CSI-A2010

Bolt-in Hardware Kit with Tether Strap
(for vehicle manufactured before 2002)

CSI-A2031

2" Extensor Thrust Wedge

CSI-A2050

Buckle Guard

ZZ-A2075

Retainer Clip Guard

CSZ-A2085