

## Folding Steel Commode

HCPCS: E0163



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Easily opens and folds
- Folds flat for convenient storage and transportation (Figure A)
- Durable plastic snap-on seat and lid (Figure B)
- Assembles without pushing buttons (Figure B - D)
- Grey powder-coated steel-welded construction increases strength and durability
- Plastic arms for added comfort
- Comes complete with 8 qt commode bucket with carry handle, cover and splash shield

PRODUCTS	Description
11148N-4	Commode,Folding,with DeepSeat,4/cs

  

SPECIFICATIONS	Description
Carton Shipping Weight	59 lbs.
Height	15.5" - 21.75"
Seat Dimensions	13.25" (W) x 16.5" (D)
Warranty	Limited Lifetime
Weight	13 lbs.
Weight Capacity	350 lbs
Width	18" between arms   22.25" Outside