



Etac Clean Height Adjustable

In addition to all the features of the basic model, this variant allows you to adjust the height of the seat.

Adjust the height for a perfect fit

The possibility of adjusting the seat height means that Etac Clean Height Adjustable can be tailor-made to suit the user's needs – for example to make transfers easier or to fit perfectly over the toilet. Adjustment is simple and no tools are required.

Etac Clean Height Adjustable has the same basic benefits as the other members of the Clean family – but is also height-adjustable up to 60 cm, and has a greater range of foot supports and calf supports.

Accessories

See our wide range of accessories "Accessories for all Clean models" on page 24–25.



Etac Clean Height Adjustable

4 lockable wheels, white

Item no.

80229208

Weight

14.2 kg

Maximum user weight

130 kg

Materials

Chassis: Powder coated steel tubing
Seat/Back support: Polypropylene
Arm supports: Polypropylene
Foot support: Polypropylene
Castors with ball, bearings: Polyamide/TPE/Stainless steel tubing
Castor height adjustment tubes: Stainless steel



Maintenance



Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max 85°C / 185°F if so required.

