

Sling Selection

Created from the online tool



General Purpose Clip Sling with Padded Legs

With padded legs for extra comfort, the general purpose sling is used for a range of everyday patient handling routines, such as bed and chair related transfers.

Article number: MAA4000M

SWL: 272 kg/600 lbs

Sizes: XS, S, M, L, LL, XL, XXL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS



General Purpose Clip Sling with Unpadded Legs

The unpadded sling is suitable for everyday patient handling routines, such as bed and chair related transfers. With no padding in the leg area, it is easier to apply to residents/patients with skeletal or spinal conditions.

Article number: MAA2000M

SWL: 272 kg/600 lbs

Sizes: XS, S, M, L, LL, XL, XXL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS

Sling Selection

Created from the online tool



General Purpose Clip Sling with Padded and Extended Legs

The extended leg sling is suitable when a larger circumference is needed to reach around the thighs. The general purpose clip slings provide good comfort and support for a range of everyday resident/patient handling routines.

Article number: MAA4100M

SWL: 272 kg/600 lbs

Sizes: XS*,S,M,L,LL,XL,XXL*

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service



General Purpose All Day Clip Sling with Unpadded Legs

Designed with the unique ability to remain under the patient between transfers while creating minimum discomfort for the resident/patient, the breathable black silk like material makes a discreet solution in combination with wheelchairs.

Article number: MAA2040M

SWL: 190 kg/418 lbs

Sizes: XS*,S,M,L,LL,XL,XXL*

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service

Sling Selection

Created from the online tool



General Purpose Comfort Deluxe Clip Sling

The Comfort Deluxe™ sling offers a very high level of comfort in key pressure areas, such as the legs or shoulders, due to no binding and additional padding. Adjustable clips give excellent head support without the use of plastic stiffeners.

Article number: MAA4050M

SWL: 228 kg/500 lbs

Sizes: XS*, S, M, L, LL, XL, XXL*

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service



Mesh Clip Sling with Padded Legs

With padded legs for extra comfort, this sling is suitable for bathing and showering as the mesh material allows water to pass easily through the sling.

Article number: MAA4061M

SWL: 272 kg/600 lbs

Sizes: XS, S, M, L, LL, XL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS

Sling Selection

Created from the online tool



Mesh Clip Sling with Unpadded Legs

Mesh material allows water to pass easily through the sling, making it an ideal choice for showering or bathing.

Article number: MAA4060M

SWL: 272 kg/600 lbs

Sizes: XS, S, M, L, LL, XL, XXL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS

Mesh Clip Sling with Extended Legs

The extended leg version has been designed to create a longer reach around the thighs. The mesh material makes it suitable for bathing allowing water to pass easily through the sling.

Article number: MAA4160M

SWL: 272 kg/600 lbs

Sizes: XS*, S*, M, L, LL, XL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service

Sling Selection

Created from the online tool



Mesh All Day Clip Sling with Unpadded Legs

An All-Day Mesh sling for general transfers is suitable for residents/patients who are in a wheelchair for a long period of time every day, and it also meets their showering or bathing needs. The mesh material allows water to pass easily through the sling.

Article number: MAA2070M

SWL: 190 kg/418 lbs

Sizes: S, M, L, LL, XL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS



Toilet Clip Sling

The toilet sling has been designed to ensure dignified handling routines with easy and efficient access to residents'/patients' clothing for dressing and undressing in connection with toileting routines.

Article number: MAA4031M

SWL: 272 kg/600 lbs

Sizes: XS, S, M, L, LL, XL, XXL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS

Sling Selection

Created from the online tool



General Purpose Paediatric Clip Sling

Every aspect of this general purpose pediatric sling has been designed with children in mind.

Article number: MAA4020M

SWL: 125 kg/275 lbs

Sizes: XXS, XS, S

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS



Double Amputee Clip Sling

Suitable for residents/patients with bilateral leg amputations above knee or total hip, this sling provides good comfort and support for a range of everyday resident/patient handling routines, such as bed and chair related transfers.

Article number: MAA2090M

SWL: 272 kg/600 lbs

Sizes: XS, S, M, L, LL, XL, XXL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Power DPS

Sling Selection

Created from the online tool



Double Amputee Clip Sling with Velcro® Fastener Flap

Suitable for residents/patients with leg amputations above or below knee, this sling is applied like a normal sling, without “log rolling”. The adjustable strap allows access to clothing for hygiene requirements and toileting.

Article number: MAA4040M

SWL: 272 kg/600 lbs

Sizes: XS*, S*, M, L, LL*, XL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service



Right Leg Amputee Clip Sling

Designed for residents/patients with right leg amputations above or below knee or total hip, this sling is shaped to give additional support around the amputation site.

Article number: MAA4070M

SWL: 272 kg/600 lbs

Sizes: S*, M, L, LL, XL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service

Sling Selection

Created from the online tool



Left Leg Amputee Clip Sling

Designed for residents/patients with left leg amputations above or below knee or total hip, this sling is shaped to give additional support around the amputation site.

Article number: MAA4080M

SWL: 272 kg/600 lbs

Sizes: S*, M, L, LL, XL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service



Double Amputee Mesh Clip Sling

Suitable for bathing and showering of residents/patients with bilateral leg amputations above knee or total hip. The mesh material allows water to pass easily through the sling.

Article number: MAA7090M

SWL: 272 kg/600 lbs

Sizes: XS*, S*, M, L, LL*, XL, XXL*

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service

Sling Selection

Created from the online tool



Double Amputee Mesh All Day Clip Sling

Suitable for bathing and showering of residents/patients with bilateral leg amputations above knee or total hip who are in a wheelchair for a long period of time every day.

Article number: MAA2080M

SWL: 190 kg/418 lbs

Sizes: XS*, S, M, L, LL, XL*

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service



Double Amputee Clip Sling All Day

Designed with the unique ability to remain under the patient between transfers while creating minimum discomfort for the resident/patient. Suitable for residents/patients with bilateral leg amputations above knee or total hip.

Article number: MAA2050M

SWL: 190 kg/418 lbs

Sizes: XS*, S, M, L, LL, XL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service

Sling Selection

Created from the online tool



Double Amputee Toilet Clip Sling

Designed with a commode hole to assist with toileting and hygiene, this sling is suitable for general transfers of residents/patients with bilateral leg amputations above knee or total hip.

Article number: MAA2091M

SWL: 272 kg/600 lbs

Sizes: S*, M, L, LL*, XL, XXL*

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Power DPS

*Available through the ArjoHuntleigh Tailor Made Sling Service



General Purpose Loop Sling without Head Support

Suitable for residents/patients who have trunk and head stability and in cases where the resident's/patient's arms, due to a medical condition or personal preference, need to be positioned outside the sling.

Article number: MLA3000

SWL: 272 kg/600 lbs

Sizes: S, M, L, XL

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point

Sling Selection

Created from the online tool



General Purpose Loop Sling with Head Support

Suitable for everyday resident/patient handling routines. The padded leg area of the sling provides high comfort helping to prevent the material roping and creasing under the resident's/patient's legs.

Article number: MLA2000

SWL: 272 kg/600 lbs

Sizes: XS*, S, M, L, XL, XXL*

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point

*Available through the ArjoHuntleigh Tailor Made Sling Service



Toilet Loop Sling without Head Support

The toilet sling without head support has been designed to ensure dignified handling routines with easy and efficient access to residents'/patients' clothing for dressing and undressing in connection with toileting routines.

Article number: MLA4531

SWL: 272 kg/600 lbs

Sizes: S, M, L, XL

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point

Sling Selection

Created from the online tool



Toilet Loop Sling with Head Support

The Toilet Sling with head support is a purpose-built design especially for toileting, dressing and undressing resident/patients.

Article number: MLA4031

SWL: 272 kg/600 lbs

Sizes: S, M, L, XL

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point



Mesh Loop Sling with Head Support

The mesh sling with head support is suitable for bathing and showering as the mesh material allows water to pass easily through the sling.

Article number: MLA4060

SWL: 272 kg/600 lbs

Sizes: XS*, S, M, L, XL, XXL*

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point

*Available through the ArjoHuntleigh Tailor Made Sling Service

Sling Selection

Created from the online tool



General Purpose All Day Loop Sling with Head Support

An all day sling is suitable for residents/patients who are in a wheelchair for a long period of time every day.

Article number: MLAAS2000

SWL: 190 kg/418 lbs

Sizes: XS*, S, M, L, XL, XXL*

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point

*Available through the ArjoHuntleigh Tailor Made Sling Service



Hammock Amputee Loop Sling

The Hammock Amputee Sling provides good comfort and support for a range of everyday resident/patient handling routines, such as bed and chair related transfers.

Article number: MLA7000

SWL: 272 kg/600 lbs

Sizes: S, M, L, XL

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point

Sling Selection

Created from the online tool



Walking Jacket with loop attachment

The Walking Jacket with loop attachment is used for standing and walking practice with a patient lift.

Article number: MAA5000

SWL: 190 kg/418 lbs

Sizes: XXS*, XS*, S, M, L, XL, XXL*

Compatible lifts and spreader bars:

Maxi Move: 2 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point

*Available through the ArjoHuntleigh Tailor Made Sling Service



Bariatric Deluxe Sling

The Bariatric Deluxe Sling provides reinforced material and stitching in the most stressed areas of the sling, including the shoulders, legs and spine.

Article number: MAA8000

SWL: 454 kg/1000 lbs

Sizes: M, L, XL, XXL

Compatible lifts and spreader bars:

Tenor: 4 point

Maxi Sky 1000: 4 point

Maxi Sky 2 Plus: 4 point (Twin Cassette)

Sling Selection

Created from the online tool



Bariatric Basic Sling

The Bariatric Basic Sling has the same shape as the Deluxe Sling, but a more economic design. It provides reinforced material and stitching in the most stressed areas of the sling, including the shoulders, legs and spine.

Article number: MAA8010

SWL: 454 kg/1000 lbs

Sizes: M, L, XL, XXL

Compatible lifts and spreader bars:

Tenor: 4 point

Maxi Sky 1000: 4 point

Maxi Sky 2 Plus: 4 point (Twin Cassette)



Bariatric Hammock Sling Div Legs

The Bariatric Hammock Sling with Divided Legs is suitable for patients needing greater stomach support. With additional support around the leg area – each leg is supported individually to prevent abduction or adduction and provides more space around the abdomen area.

Article number: MAA8020

SWL: 454 kg/1000 lbs

Sizes: M, L, XL, XXL

Compatible lifts and spreader bars:

Tenor: 4 point

Maxi Sky 1000: 4 point

Maxi Sky 2 Plus: 4 point (Twin Cassette)

Sling Selection

Created from the online tool



Bariatric Hammock Sling

The Bariatric Hammock Sling is suitable for transfers when patient is semi reclined and may also be used with amputee patients. May be suitable for patients who cannot tolerate leg straps positioned between their knees.

Article number: MAA8030

SWL: 454 kg/1000 lbs

Sizes: M, L, XL, XXL

Compatible lifts and spreader bars:

Tenor: 4 point

Maxi Sky 1000: 4 point

Maxi Sky 2 Plus: 4 point (Twin Cassette)



Bariatric Disposable Repositioning Sling

The patient specific Bariatric Repositioning Sling assists the caregiver in bed repositioning and transferring heavy patients in a supine position between beds/stretchers. The sling dimensions fit with the bariatric sized beds.

Article number: VIG22004-L1

SWL: 454 kg/1000 lbs

Sizes: One size 145 x 210 cm, 57 x 82 "

Compatible lifts and spreader bars:

Maxi Sky 1000: 4 point

Maxi Sky 2 Plus: 4 point (Twin Cassette)

Sling Selection

Created from the online tool



Bariatric Repositioning Sling

Assisting the caregiver in bed repositioning and transferring heavy patients in a supine position between beds/stretchers. The sling dimensions fit with the bariatric sized beds.

Article number: VIG11004-L1

SWL: 454 kg/1000 lbs

Sizes: One size 145 x 210 cm, 57 x 82 "

Compatible lifts and spreader bars:

Maxi Sky 1000: 4 point

Maxi Sky 2 Plus: 4 point (Twin Cassette)



Bariatric Limb Sling

Designed with a wider undisturbed surface and higher safe working load than the standard limb slings.

Article number: VIG11002-L1

SWL: 120 kg/264 lbs (454 kg/1000 lbs full patient weight)

Sizes: One size 180 x 30 cm, 70 x 12"

Compatible lifts and spreader bars:

Tenor: 4 point

Maxi Sky 600: 2, point

Maxi Sky 1000: 2, point 4 point

Maxi Sky 2: 2 point, 4 point

Sling Selection

Created from the online tool

Maxi Transfer Sheet

The Maxi Transfer Sheet is a unique and innovative combination of bed linen and sling, setting a new standard for lateral and in-bed positioning solutions.

Article number: MAA7000

SWL: 272 kg/600 lbs

Sizes: Regular, Wide

Compatible lifts and spreader bars:

Maxi Move: 2 point, Stretcher Frame

Maxi Sky 600: 2 point, Stretcher Frame

Maxi Sky 2: 2 point, 4point, Stretcher Frame



Soft Stretcher Mesh Sling with Commode Hole

With mesh fabric and commode opening, the Soft Stretcher Sling is suitable for hygiene of residents/patients who require full body support in a supine position during transfers.

Article number: MAA1152

SWL: 272 kg/600 lbs

Sizes: L, XL

Compatible lifts and spreader bars:

Maxi Move: Stretcher Frame

Maxi Sky 600: Stretcher Frame

Maxi Sky 2: Stretcher Frame

Sling Selection

Created from the online tool



Soft Stretcher Sling with Commode Hole

Including a commode opening, this sling is suitable for residents/patients who require full body support in a supine position during transfers.

Article number: MAA1150

SWL: 272 kg/600 lbs

Sizes: L, XL

Compatible lifts and spreader bars:

Maxi Move: Stretcher Frame

Maxi Sky 600: Stretcher Frame

Maxi Sky 2: Stretcher Frame



Soft Stretcher Sling without Hole

Suitable for residents/patients who require full body support in a supine position during transfers.

Article number: MAA1151

Sizes: L, XL

Compatible lifts and spreader bars:

Maxi Move: Stretcher Frame

Maxi Sky 600: Stretcher Frame

Maxi Sky 2: Stretcher Frame

Sling Selection

Created from the online tool



Repositioning Sling

Repositioning Sling for regular repositioning in bed. It is mainly used for moving the resident/patient up the bed or turning them.

Article number: MAA6000

SWL: 272 kg/600 lbs

Sizes: One size 204 x 104 cm, 80 x 41"

Compatible lifts and spreader bars:

Maxi Move: 2 point

Maxi Twin: 2 point

Maxi 500: 2 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point



Limb Loop Sling

Designed to assist the caregiver with lifting, supporting and positioning the limb for easier access when performing activities such as dressing, foot care, wound care- and venous wrapping.

Article number: VIG11001-L1

SWL: 70 kg/154 lbs (272 kg/600 lbs full body weight)

Size: One size: 120 x 14 cm, 47 x 5,5 "

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Maxi Lite: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point



Limb Clip Sling

Compatible with ArjoHuntleigh Clip Lift Systems, this sling is designed to assist the caregiver with lifting, supporting and positioning the limb for easier access when performing activities such as dressing, foot care, wound care and venous wrapping.

Article number: VIG11000-L1

SWL: 70 kg/154 lbs (272 kg/600 lbs full body weight)

Size: One size: 120 x 14 cm, 47 x 5,5 "

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS



Turning Sling

This sling is designed to be applied without having to log roll or move the patient and assists the caregiver with turning/rolling a patient in a supine position in bed. The narrow surface area of the sling allows access to patient's skin for inspection, treatment and personal hygiene.

Article number: VIG11003-L1

SWL: 100 kg/220 lbs (272 kg/600 lbs full patient weight)

Size: One size: 80 x 45 x 14 cm, 31,5 x 17,7 x 5,5 "

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Maxi Lite: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point

Sling Selection

Created from the online tool



Sara 3000 Standing Sling with Fixed Support Belt

Sara® 3000 Standing Sling for all general transfers and toileting activities with a fixed belt.

Article number: TSS.50X

SWL: 200 kg/440 lbs

Sizes: S, M, L, XL, XXL

Compatible lifts:

Sara 3000
Sara Lite
Sara Steady



Sara 3000 Standing Sling with Separate Support Belt and 1 Clip Closing

The separate support belt is narrower than the fixed belt and attached with one belt clip instead of two. Suitable for all general transfers and toileting activities.

Article number: TSS.50XS-1C

SWL: 200 kg/440 lbs

Sizes: S*, M, L, XL, XXL*

Compatible lifts:

Sara 3000
Sara Lite
Sara Steady

*Available through the ArjoHuntleigh Tailor Made Sling Service

Sling Selection

Created from the online tool



Sara 3000 Standing Sling with Separate Support Belt and Velcro® Fastener

The belt is attached with Velcro® fasteners instead of belt clips and narrower than the fixed belt. Suitable for all general transfers and toileting activities.

Article number: TSS.50XSV

SWL: 200 kg/440 lbs

Sizes: S*, M, L, XL*, XXL*

Compatible lifts:

Sara 3000
Sara Lite
Sara Steady

*Available through the ArjoHuntleigh Tailor Made Sling Service



Sara 3000 Transfer Sling

Sara® 3000 Transfer Sling for supported transfers in a seated position.

Article number: TSS.511

SWL: 200 kg/440 lbs

Sizes: One size

Compatible lifts:

Sara 3000

Sling Selection

Created from the online tool

ARJOHUNTLEIGH
GETINGE GROUP



Sara Plus Standing Sling

Sara® Plus Standing Sling for all general sit to stand, transfer and toileting activities.

Article number: KKA5120

SWL: 190 kg/418 lbs

Sizes: S, M, L, XL

Compatible lifts:
Sara Plus



Sara Plus EPS Sling

Extra Postural Support (EPS) to help patients achieve full extension when standing. Suitable for all general sit to stand, transfer and toileting activities.

Article number: KKA5370

SWL: 190 kg/418 lbs

Sizes: S, M, L, XL

Compatible lifts:
Sara Plus

Sling Selection

Created from the online tool



Sara Plus Wipe Down Sling

Designed with a wipeable material this sling can be wiped clean and disinfected between patients. Suitable for all general sit to stand, transfer and toileting activities.

Article number: KKA5420

SWL: 190 kg/418 lbs

Sizes: S, M, L, XL

Compatible lifts:
Sara Plus



Sara Plus Walking Sling

Suitable when Sara Plus is used for standing, stepping and walking practice rather than the traditional use of a standing and raising aid, this sling gives safe support to patients during standing, stepping and walking practice.

Article number: KKA5130

SWL: 190 kg/418 lbs

Sizes: S, M, L, XL

Compatible lifts:
Sara Plus

Sling Selection

Created from the online tool



General Purpose Clip Flites

General Purpose Clip Flites® is designed to be used for one patient during a hospital stay and assists the caregiver with everyday patient handling routines, such as bed and chair related transfers.

Article number: MFA1000M

SWL: 272 kg/600 lbs

Sizes: S, M, L, LL, XL, XXL

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS



General Purpose Loop Flites

General Purpose Loop Flites® combines infection control and safe patient handling routines. It is suitable for everyday patient handling routines, such as bed and chair related transfers.

Article number: MFA2000

SWL: 272 kg/600 lbs

Sizes: S, M, L, XL

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point

Sling Selection

Created from the online tool



Sara 3000 Flites

The patient specific Sara® 3000 Flites combines infection control and safe patient handling routines. It's suitable for all general transfers and toileting activities in a standing position. The fixed belt is attached with two belt clips.

Article number: MFA3000

SWL: 200 kg/440 lbs

Sizes: One size

Compatible lifts:

Sara 3000

Sara Lite

Sara Steady



Sara Plus Flites

Sara® Plus Flites for all general sit to stand, transfer and toileting activities with the Sara Plus standing and raising aid combining infection control and safe patient handling routines.

Article number: MFA4000

SWL: 190 kg/418 lbs

Sizes: One size

Compatible lifts:

Sara Plus

Sling Selection

Created from the online tool



Stretcher Flites

Patient specific Stretcher Flites® for patients who require full body support in a supine position during transfers. This full body one piece sling has been designed for use with a stretcher frame.

Article number: MFA5000

SWL: 272 kg/600 lbs

Sizes: L

Compatible lifts and spreader bars:

Maxi Move: Stretcher Frame

Maxi Sky 600: Stretcher Frame

Maxi Sky 2: Stretcher Frame



Limb Clip Flites

Designed for single patient use the Limb Clip Flites assists the caregiver with lifting, supporting and positioning the limb for easier access when performing activities such as dressing, foot care, wound care- and venous wrapping. Compatible with ArjoHuntleigh Clip Lift Systems.

Article number: VIG22000-L1

Size: One size: 120 x 14 cm, 47 x 5.5"

SWL: 70 kg/154 lbs (272 kg/600 lbs full patient weight)

Compatible lifts and spreader bars:

Maxi Move: Manual DPS, Power DPS

Maxi Twin: Manual DPS, Power DPS

Maxi Twin Compact: Power DPS

Maxi Lite: Flat DPS

Maxi 500: Manual DPS

Maxi Sky 2/600: Manual Flat DPS, Power DPS



Limb Loop Flites

Designed for single patient use the Limb Loop Flites assists the caregiver with lifting, supporting and positioning the limb for easier access when performing activities such as dressing, foot care, wound care- and venous wrapping.

Article number: VIG22001-L1

SWL: 70 kg/154 lbs (272 kg/600 lbs full body weight)

Size: One size: 120 x 14 cm, 47 x 5,5 "

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Maxi Lite: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point, 4 point



Turning Flites

The patient specific disposable Turning Flites assists the caregiver with turning a patient in a supine position in bed. The sling is designed to be applied without having to log roll or move the patient.

Article number: VIG22003-L1

SWL: 100 kg/220 lbs (272 kg/600 lbs full patient weight)

Size: One size: 80 x 45 x 14 cm, 31,5 x 17,7 x 5,5 "

Compatible lifts and spreader bars:

Maxi Move: 2 point, 4 point

Maxi Twin: 2 point, 4 point

Maxi 500: 2 point

Maxi Lite: 2 point

Minstrel: 4 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point

Maxi Sky 2: 2 point

Sling Selection

Created from the online tool



Disposable Repositioning Sling

The patient specific Disposable Repositioning Sling assists the caregiver in bed repositioning and moving the patient up in bed . Regular repositioning, together with pressure area care, provides efficient pressure management and thereby helps to prevent pressure ulcers.

Article number: AHD001

SWL: 272 kg/600 lbs

Sizes: One size 208 cm x 119 cm, 82 x 47 "

Compatible lifts and spreader bars:

Maxi Move: 2 point

Maxi Twin: 2 point

Maxi 500: 2 point

Maxi Sky 440: Fixed 2 point

Maxi Sky 600: 2 point