



# R82 Cricket

Compact and Lightweight  
Easily Foldable Wheelchair

R82<sup>®</sup> 



# R82 Cricket

The R82 Cricket is a lightweight, compact and easy-to-fold wheelchair developed with aesthetics, quality and comfort in mind. This fun push wheelchair provides support and stability for the child, along with excellent functionality for the caregivers. With its effortless maneuverability, the R82 Cricket is the perfect chair for outings and adventures! This 20° fixed tilt chair comes in two sizes and adapts to the child's growth, maximizing its lifespan. Children will delight in the choice of colors and caregivers will appreciate the ease of use and flexibility of the R82 Cricket. The R82 Cricket is approved for safe and portable transportation and adheres to ISO 7176-19 standards.





## Convenience and Ease of Use

- **Tool-less Back Recline (110° - 125°)**  
The backrest angle can be changed by 15° with a simple upholstery release
- **Lightweight and Easy to Fold**  
This lightweight wheelchair is compact when folded and is self-standing, requiring minimum storage space
- **Growth Capabilities**  
As the user grows, Cricket's seat depth can grow 4" and height of footplate can be adjusted according to the child's size



Isabella in an R82 Cricket walking with her family.

## Accessories



### Armrests

Size 1: H 6" x L 8" x W 2"  
Size 2: H 7.5" x L 10.5" x W 2.5"



### Back and Seat Contour Wedges

Size 1: W 4.25" x L 5.75" x H 1"  
Size 2: W 4.75" x L 5.75" x H 1"



### Backrest Extension

Size 1: H 8.5" x W 10"  
Size 2: H 8.5" x W 12"




### Seat Width Reduction Wedges

Size 1: H 5.75" x L 6" x W 1"  
Size 2: H 5.75" x L 6" x W 2"



### Calf Panel

Size 1: L 14.5" x W 4.5"  
Size 2: L 17" x W 4.5"



### Underseat Storage Basket

Size 1: L 14" x W 15"  
Size 2: L 14" x W 17.5"



### Canopy

Size 1: H 16.5" x W 14"  
Size 2: H 18.5" x W 16"



### Sculpted Headrest

H 7" x W 11.25" x D 7"



### Headrest Cushion

H 5.75" x W 11" x D 3.5"




### Tray

L 16" x D 13.5"



### Heel Straps

H 3" x L 7"




### Grip Bar

In frame



### Grip Bar

In armrests





















### Rain Cover

Please visit [R82.com](http://R82.com) to see additional options, including the vests and belts that can be used with the R82 Cricket.


# R82 Cricket Specifications

inches

	1	2
 Seat width	13½"	15¾"
 Seat depth	8½–12¼"	12¼–17"
 Back height	19¾"	21¼"
 Seat to floor height	19¾–21"	21¼–22¼"
 Seat to footplate distance	7–14¼"	8¾–14½"
 Leg to seat surface angle	360°	360°
 Backrest width	12¾"	15¼"
 Backrest width, inside	10¾"	13"
 Overall width	26¾"	29¼"
 Back recline	110° - 125°	110° - 125°

	1	2
 Overall length	40¼"	41¼"
 Overall height	39¼"	41¼"
 Folded width	14"	14"
 Total weight	27.3 lbs	29.5 lbs
 Max. user weight/load	77 lbs	110 lbs
 Max. user weight/load, transportation	59.5 lbs	94.5 lbs
 Required width of right angled corridor	32"	32"
 Max. Static stability	10° sideways	10° forward 10° rearward

## Cushion Color Options



Everyone should be able to live a free and independent life  
and pursue their dreams regardless of any physical circumstances.

R82 Inc.  
12801 E. Independence Blvd.  
Matthews, NC 28106  
Tel +1 844-US-MOBILITY (844-876-6245)  
Fax +1 704-882-0751  
www.R82.com  
information@R82.com  
RB1001 0717

