



## FEATURES

- Orthopedic contoured backs in Hi-Back™ and Lo-Back™ options (Hi-Back models include harness for added support)
- Stainless steel heavy duty mounting base
- Padded back for pressure reduction
- Latex and silicone free

Shown with Optional  
Armrest and Footrest



Shown with Reducer Ring



U-Shape Mounting Base allows you to upgrade or replace the support back to the next size for a fraction of the cost



Reduced seat depth (RSD): sits 4" forward to lessen seat depth when attached to U-Shape Mounting Base

## ACCESSORIES



Armrest & Footrest



Padded Reducer Ring



Plastic Reducer Ring



Soft-Flex Splash Guard



**ORDER FORM** **FAX TO 310-305-1718 OR SUBMIT ONLINE**

Dealer: \_\_\_\_\_ Marked For: \_\_\_\_\_ Account #: \_\_\_\_\_

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Ordered By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_ E-mail: \_\_\_\_\_

Ship to address: \_\_\_\_\_

<b>LO-BACK BASE UNIT</b>	<b>TCL-4210S*</b>	<b>TCL-4310M</b>	<b>TCL-4410L</b>	
Product Size	Small	Medium	Large	
Inside Back Width at Top	7"	11"	13"	
Inside Back Width at Bottom	7"	11"	13"	
Back Height Range	7" - 14"	10" - 14"	10" - 17"	
Product Weight	5.5 lb	6.5 lb	7.75 lb	
<b>HI-BACK BASE UNIT</b>	<b>TCH-4610S*</b>	<b>TCH-4710M*</b>	<b>TCH-4810M</b>	<b>TCH-4910L</b>
Product Size	Small	Medium	Medium	Large
Inside Back Width at Top	7"	8"	8"	11"
Inside Back Width at Bottom	12"	12"	12"	14"
Back Height Range	17" - 22"	19" - 24"	19" - 24"	24" - 28"
Product Weight	9.5 lb	9.75 lb	9.75 lb	12 lb

\* Reduced seat depth (RSD): sits 4" forward to lessen seat depth when attached to U-Shape Mounting Base

**ACCESSORIES**

Attached Armrest	TZ-A4921		
Footrest (must use with armrest)	TZ-A4922		
Plastic Reducer Ring (Hole: 6"W x 7"L)	TZ-A4950		
Padded Reducer Ring with Abductor	TZ-A4965XS	Hole: 5"W x 6"L	
	TZ-A4960S	Hole: 7"W x 8"L	
Soft-Flex Splash Guard (5"W x 5.5"L)	TZ-A4975		
U-Shape Mounting Base (18.5"W x 3"D x 10"H)	TZ-R4001		