



# ***MallMASTER SPRINT S237***



- 3 Wheel Scooter
- Lightweight & Portable
- Deluxe Swivel Sliding Seat
- Flat Free Tires
- 22 ah Detachable Battery Pack
- Colours: Red, Blue and White

		<b>S237</b>
Dimensions L x W x H	cm (in)	102 (40) x 54 (21) x 82 (32)
Total Weight w/Batteries	kg (lb.)	47 (104)
Battery	_V_Ah x Pcs	12V / 22 AH x 2
Charger	Off Board	1.8 amp
Tires		Front: 8" / Rear: 9"
Motor	HP / DC24V	.6 / 480W (Max)
Max. Speed	km/hr. (mph)	7.2 (4.5)
Driving Distance*	km (mile)	17 (10.6)
Weight Capacity	kg (lb.)	136 (300)
Suspension		No
Turning Radius	cm (in)	86 (34)

Ver 1.2

Specifications subject to change without notice

Product may not be exactly as shown.

\*Maximum driving distance is based on a new fully charged battery depending on weight, terrain and temperature.

Discovering your world again—in style

