

Bariatric Transfer Bench

HCPCS: E0248



FIGURE A



FIGURE B

- Tool-free assembly of back, legs and arm
- A-frame provides additional stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- All transfer benches are reversible to accommodate any bathroom
- Euro-style clip ensures a secure fit (Figure A)
- No exposed hardware to injure patient
- Large, locking suction cups for added safety (Figure B)
- Weight Capacity 500 lbs

PRODUCTS

Description

12025KD-1	Transfer Bench, Heavy Duty, KD, 1/cs
-----------	--------------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	14 lbs. 27 lbs.
Outside Legs	32.5" (W) x 24" (D)
Seat Dimensions	26.5" (W) x 16.5" (D)
Seat Height	18.5" - 23.5"
Warranty	Limited Lifetime
Weight	14 lbs.
Weight Capacity	500 lbs

PARTS

Description

12011SS	SuctionCups,Sm.12011,1pr/bx, BX
12011SL	SuctionCups,Lg.12011,1pr/bx, BX
12021C	Horseshoe Clip,Ht Adjust, 12021, 1/ea, EA
12025KDC	Back Clip,12025KD,1ea, EA
12025KDB	Back Assembly,12025KD,1set/pkg, PKG
12025KDLGGS	Set of 4 Legs,12025KD,1set/pkg, PKG
12025KD9LB	Back, 12025KD, 9L Ser#, 1/ea, EA
12025KD9LS	Seat, 12025KD, 9L Ser#, 1/ea, EA