

Fusion

APT Series

Manufactured by



Alternating Pressure

The Fusion APT system features transfer edges, a quilted cover, and quiet operation. This results in a more restful and restorative sleep, as well as safer transfers onto an air surface. The mattress has a 5" air bladder that provides therapeutic alternating pressure with true offloading. The Fusion APT has an added feature that allows the mattress to stay inflated in a power outage for up to 24 hours.

Advanced Pump Options

With 8 adjustable levels of comfort, you get the settings that work best for you and your therapy requirements. The Fusion pump also has a 2-in-1 alternating 10-minute cycle to ensure maximum protection. For peace of mind, a low-pressure alarm is triggered when the air pressure of the mattress is lower than the default setup.

Safer Transfers

The 3" wide transfer edges on the Fusion APT enables safer transfers on and off the surface, and reduces the risk of falls and bed entrapment.



Specifications

Cover Type	Top: Quilted Recovery 5 cover Base: Polyester laminated with PVC
Size	80" x 36" x 6" (203 cm x 91.4 cm x 15.2 cm)
Product Weight	20.6 lb (9.3 kg)
Patient Weight Limit	350 lb (158.8 kg)
Warranty	36-Month Warranty