

Pinniped Pediatric Commode




FIGURE A



FIGURE B

- All-in-one welded steel frame available in red and blue
- Plastic seat and contoured back
- Durable snap-on toilet seat
- Height adjustable seat allows child to firmly place their feet on the floor for comfort, safety and security
- Adjustable hook & loop chest strap for additional support
- Molded plastic armrests provide additional comfort
- Rubber non-skid tips ensure that the legs will not mark the floor
- Complete with 7.5 qt removable commode bucket with metal handle and cover

PRODUCTS

Description

PC 1000 BL	Pinniped Pediatric Commode, Blue, 1/ea
PC 1000 R	Pinniped Pediatric Commode, Red, 1/ea

SPECIFICATIONS

Description

Armrest Height	15.5" - 19.5"
Armrest Pad	2" (W) x 9" (D)
Carton Shipping Weight	13 lbs
Overall Width	17"
Seat Depth	15"
Seat Height	10.5" - 14.5"
Seat Width	13"
Toilet Opening	6.5" (W) x 8" (D)
Warranty	Limited Lifetime
Weight	9.5 lbs
Weight Capacity	250 lbs