

3-Wheel

Model Number	S25
Weight Capacity	113.4 kg
Maximum Speed ¹	Up to 6 km/h
Ground Clearance ³	1.3"
Turning Radius ³	33.5"
Overall Length ³	38"
Overall Width ³	19" - 20"
Seat-to-Ground Height ³	13" – 15"
Seat-to-Deck Height ³	8.75" – 10.75"
Front Tires ³	6.7" solid
Rear Tires ³	6.7" solid
Suspension	N/A
Range Per Charge ^{1,2}	Up to 10.8 km/h
Total Weight Without Batteries ⁵	22.2 kg (21 kg with batteries)
Weight of Heaviest Piece ³	14.6 kg (Base only)
Standard Seating	Type: Compact Weight: 4.4 kg Material: Molded plastic Dimensions: 16" x 16"
Drive System	36 volt brushless hub motor
Dual Braking System	Manual, dual disc
Battery Requirements ^{4,6}	Size: (1) 36 volt, 5.2 AH (lithium ion) Weight: 2 kg
Battery Charger	36 volt lithium, off-board
Warranty	Frame: lifetime limited Drivetrain: 1-year limited Electronics: 1-year limited Battery: 6-months

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

⁵ Includes standard seat weight.

⁶ AGM or gel cell type required.



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your iCruise and its lithium ion battery are permitted.

PRIDE
MOBILITY
LIVE YOUR BEST®

Beamsville, ON • 888-570-1113

pridemobility.ca |    

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

GO-GO_ICRUISE_040120

PRELIMINARY

GO-GO
ICRUISE

PORTABLE.
LIGHTWEIGHT.
iCRUISE.



LIVE YOUR BEST®

pridemobility.ca |    

LIGHTWEIGHT DESIGN, EFFORTLESS TRAVEL



Digital LCD display




Manual hand brake



Supplemental foot brake



Twist grip throttle



Removable seat with foldable tiller for effortless portability

Journey far and wide with the **iCruise!** This ultra-lightweight scooter is super portable and perfect for traveling. Equipped with both a manual hand brake and supplemental foot brake, the iCruise delivers a safe and comfortable ride!

- Portable and easy to travel
- Supplemental foot brake
- Manual hand brake
- Digital LCD display
- Foldable tiller with adjustment knob
- Airline-safe lithium ion battery pack

Available Colours:



Lime

Mango

Blackberry



Arctic Ice

Raspberry

Options

- Additional 5.2 AH battery