



New Nimbo 2G Walker with Color Options!

The already popular Nimbos are now available in three colors to allow children using these walkers to select a color that matches their own personal style. The improved hand grip design offers smoother texture for increased comfort.



SAME GREAT FEATURES

- Height adjustable aluminum frame
- Easily folds for storage and transport
- Soft rubber wheels adhere to any surface and prevent the Nimbo from sliding backwards when using one-directional setting
- One-directional 5" rear wheels with aluminum ratchet and pin mechanism for increased durability
- One-directional override bracket to allow for two-directional movement
- Easily switch from swivel to non-swivel front wheel settings

FRAME



Knight Blue



Castle Red



Wizard Purple



Emperor Black
(Available in XLarge Only)

ACCESSORIES



Soft Seat Harness



Pelvic Stabilizer



Forearm Platforms



Fold-down Seat



Bumper Guard



Dial-A-Speed Tab

ORDER FORM

FAX TO 905-787-8550 OR SUBMIT ONLINE













Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims. Referenced HCPCS codes are also not guarantee of payment.

BASE		Extra Small	Small	Medium	Large	Extra Large
Castle Red		KA1200-2GCR	KA2200-2GCR	KA3200-2GCR	KA4200-2GCR	N/A
Knight Blue		KA1200-2GKB	KA2200-2GKB	KA3200-2GKB	KA4200-2GKB	N/A
Wizard Purple		KA1200-2GWP	KA2200-2GWP	KA3200-2GWP	KA4200-2GWP	N/A
Emperor Black		N/A	N/A	N/A	N/A	KA5200-2GEB
OPTIONAL SEAT		KA1285-2GCR	KA2285-2GCR	KA3285-2GCR	KA4285-2GCR	N/A
		KA1285-2GCR	KA2285-2GCR	KA3285-2GCR	KA4285-2GCR	N/A
		KA1285-2GKB	KA2285-2GKB	KA3285-2GKB	KA4285-2GKB	N/A
		KA1285-2GWP	KA2285-2GWP	KA3285-2GWP	KA4285-2GWP	N/A
		N/A	N/A	N/A	N/A	KA5285-2GEB
BASE WITH SEAT		KA1200S-2GCR	KA2200S-2GCR	KA3200S-2GCR	KA4200S-2GCR	N/A
		KA1200S-2GCR	KA2200S-2GCR	KA3200S-2GCR	KA4200S-2GCR	N/A
		KA1200S-2GKB	KA2200S-2GKB	KA3200S-2GKB	KA4200S-2GKB	N/A
		KA1200S-2GWP	KA2200S-2GWP	KA3200S-2GWP	KA4200S-2GWP	N/A
		N/A	N/A	N/A	N/A	KA5200S-2GEB
Weight Capacity		75 lb	85 lb	150 lb	190 lb	250 lb
Width Inside Handlebar		13.5"	13.5"	15"	16.5"	17.75"
Handle Height		15.5" - 19.5"	19" - 25"	23" - 30.5"	28" - 35"	36" - 41"
Overall Width (open at base)		24"	24"	26"	26.5"	29"
Overall Depth (open at base)		20"	27"	31"	34"	37"
Product Weight		9 lb	10 lb	13 lb	14 lb	16 lb
Seat to Floor Height (S models)		9" - 13"	11.5" - 17"	12" - 19.5"	14.5" - 22"	21.5" - 27.5"
HCPCS		E1399	E1399	E1399	E1399	E1399

ACCESSORIES

Forearm Platforms*	KA 1035 FPS	Length: 5", Height from handle bar: 3.5" - 4.5"
	KA 1035 FPL	Length: 8", Height from handle bar: 5" - 8"
Pelvic Stabilizer	KA 1050	Inside width range: 6" - 13"
Soft Seat Harness*	CE 1070S	Can be used on all sizes of Nimbo
	CE 1070L	Can be used on all sizes of Nimbo
Dial-A-Speed Tabs	KA 1400	
Bumper Guards	KA 8000	
9" Non-Swivel Front Wheels	KA 8100	For use with Nimbo size S, M, L, XL

*Both S & L sizes can be used on all sizes of Nimbo