

# XPRESSO Series



## SPECIFICATIONS

### MINI

18 lbs\*  
18" Seat Width  
18" Seat Height  
22.5" Width x 26" Length  
30.5" - 33.5" Handle Height  
350 lbs Weight Capacity

### REGULAR

18 lbs\*  
18" Seat Width  
21" Seat Height  
22.5" Width x 26" Length  
30.5" - 34.5" Handle Height  
350 lbs Weight Capacity

### TALL

18 lbs\*  
18" Seat Width  
24" Seat Height  
22.5" Width x 26" Length  
33.5" - 40" Handle Height  
350 lbs Weight Capacity

\*Net Weight of walker without accessories, backstrap and folding basket

THE XPRESSO, ALONG WITH ALL EVOLUTION WALKERS, WAS DESIGNED TO BE THE STRONGEST AND MOST STABLE OF ALL THE CENTER FOLDING WALKERS ON THE MARKET TODAY.

With cable-free brakes, a center folding seat, durable folding basket, the Xpresso folds easily and compact.

## FEATURES

**INNOVATIVE DESIGN** Roman arch provides support & stability

**SEAT** Center folding durable waterproof polyurethane cushion

**BRAKE SYSTEM** Cable-free / Easy to adjust

**COLLAPSIBLE BASKET** Foldable & removable

**WHEELS** 8" long lasting polyurethane tires, suitable for indoor & outdoor use

**COLORS** Ruby Red, Azure Blue, Ebony Black



## OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Deluxe Back Strap Pad
- Oxygen Tank Holder
- Wire Basket
- Eyeglass Holder
- Phone Holder