




Anchor



Bittersweet



Oxford



Port



Sandstorm

Shown in Anchor

## Comforter PR531



*Comfort Zone*  
by GOLDEN

The Comforter Series from Golden is all about finding the perfect size for the perfect fit. Available in seven sizes to comfortably fit people with weight capacities up to 500 lbs. These power lift recliners feature a new back style with a tasteful seam design, our signature seat for the utmost in comfort and support, a new arm style with additional plush padding for more comfort and value, and extended arms to provide better stability while entering or exiting the chair. Enjoy luxurious new chenille fabrics to match any décor. Find your perfect fit!

### Standard Features:

- Signature seat
- Luxurious chenille fabrics
- Upgraded stain defense
- USB charging port



# Comforter Series

## PR531



| Model Number                                                   | PR531-PSA    | PR531-MED    | PR531-LAR    | PR531-TAL    | PR531-S23    | PR531-M26    | PR531-T28    |
|----------------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|
| Size                                                           | Petite/Small | Medium       | Large        | Tall         | Small/Wide   | Medium/Wide  | Tall/Wide    |
| Weight Capacity                                                | 375 LB       | 375 LB       | 375 LB       | 375 LB       | 375 LB       | 500 LB       | 500 LB       |
| Positioning                                                    | 3            | 3            | 3            | 3            | 3            | 3            | 3            |
| Comfort Zone(s)                                                | 1            | 1            | 1            | 1            | 1            | 1            | 1            |
| Twilight Technology                                            | NO           | NO           | NO           | NO           | NO           | NO           | NO           |
| Chaise Seat                                                    | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| Pocketed Coil Seat                                             | YES          | YES          | YES          | YES          | NO           | NO           | NO           |
| Back Type                                                      | Seam         | Seam         | Seam         | Seam         | Seam         | Seam         | Seam         |
| Removable Back Type                                            | Bolt         | Bolt         | Bolt         | Bolt         | Bolt         | Bolt         | Bolt         |
| Head Rest Cover                                                | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| Arm Rest Covers                                                | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| Overall Width                                                  | 31.5"        | 32"          | 34"          | 34.5"        | 34"          | 38"          | 40"          |
| Overall Height                                                 | 40"          | 41"          | 44"          | 47"          | 40.5"        | 43.5"        | 47.5"        |
| Width Between Arms                                             | 21"          | 21"          | 23"          | 23"          | 23"          | 26"          | 29"          |
| Floor to Top of Seat                                           | 19"          | 20"          | 21"          | 22.5"        | 19.5"        | 20.5"        | 22"          |
| Seat Depth                                                     | 19"          | 19"          | 20.5"        | 22"          | 20"          | 20.5"        | 22"          |
| Seat to Top of Back                                            | 27"          | 27"          | 28.5"        | 30.5"        | 27"          | 29"          | 32"          |
| Seat to Top of Arm                                             | 6"           | 5"           | 6"           | 6"           | 7"           | 7.5"         | 5.5"         |
| Distance Required from Wall Reclined                           | 16"          | 17"          | 19"          | 21"          | 16"          | 18"          | 21"          |
| Extension Reclined from Back of Seat to End of Foot Rest       | 41°          | 41"          | 42"          | 47"          | 40"          | 42"          | 46"          |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| Warranty on Electrical Parts                                   | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       |
| Pro Rated Years 4 Thru 7                                       | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| Warranty on Mechanical Labor                                   | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       | 3 Yrs.       |
| Warranty on Electrical Labor                                   | 1 Yr.        | 1 Yr.        | 1 Yr.        | 1 Yr.        | 1 Yr.        | 1 Yr.        | 1 Yr.        |
| Hand Control Routing                                           | Inside Arm   | Inside Arm   | Inside Arm   | Inside Arm   | Inside Arm   | Inside Arm   | Inside Arm   |
| Extra Pocket Option                                            | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| Heat & Massage Option                                          | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| Power Pillow Option                                            | NO           | NO           | NO           | NO           | NO           | NO           | NO           |
| Power Adjustable Lumbar Option                                 | NO           | NO           | NO           | NO           | NO           | NO           | NO           |
| AutoFlex Option                                                | NO           | NO           | NO           | NO           | NO           | NO           | NO           |
| Foot Rest Extension Option                                     | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| Straight Lift Option                                           | YES          | YES          | YES          | YES          | YES          | YES          | YES          |
| USB Charging Port                                              | Hand Control | Hand Control | Hand Control | Hand Control | Hand Control | Hand Control | Hand Control |
| Place of Origin                                                | USA          | USA          | USA          | USA          | USA          | USA          | USA          |


Battery Backup & SmartTek Diagnostic System STANDARD on all Golden models.



*Comfort Zone*  
by GOLDEN

### ONE ZONE POSITIONING SYSTEM

Comfort Zone 1 provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.



Take a picture with your smartphone to learn more.



Phone: 800-624-6374  
Fax: 800-628-5165  
[www.goldentech.com](http://www.goldentech.com)

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with  
UL 73 and  
CSA C22.2 No. 68



[www.facebook.com/goldentechnologies](http://www.facebook.com/goldentechnologies)

SHT-COMF-121820