

Invacare's line of recliners delivers outstanding value, durability and comfort in health care chairs and recliners. Our Deluxe Mobile Recumbent Chairs (MRC) provide essential ergonomic support and comfort in a recliner. The Mobile Treatment Chair (MTC) is especially suited for clinical needs. Regardless of your need, Invacare's extensive line of recliners are ideal for a range of long-term care and clinical seating applications. They're rugged, easy to clean and perform a variety of functions safely and effectively, making them an outstanding choice for you and your residents.



Traditional Three-Position Recliner

Model no. IH6074A

The IH6074A three-position recliner offers exceptional comfort and durability for day room or resident room environments. The full recline, leg-ottoman elevated and full upright positions provide maximum flexibility for residents. A convenient push bar on the chair back and four rolling casters allow for easy movement. The two rear steering casters may be locked to secure the chair in place.

Upholstery Colors



#61 Blueridge



#68 Jade



#60 Rosewood



Full recline position.

Features

- Footrests fold completely under chair, allowing the resident to rise to a standing position more easily
- Blow-molded side panels remove quickly, with press pins, for easy cleaning
- Convenient overlap tray included. Tabs allow tray to move freely
- Articulating headrest joins resident to social field of vision
- Flame retardant upholstery – meets California Technical Bulletin 117 standard

Invacare's Uni-Recliners reduce pressure, preserve independence and improve the quality of life for the resident. When a resident requires greater positioning support than a wheelchair can provide, Invacare recliners provide the solution.


Specifications

	IH6074A	IH6065	IH6077A
Overall height	47"	46"	49"
Overall width	24.25"	27.5"	30.75"
Overall upright depth	40"	41"	38"
Seat width	19"	18.75" Adult 22.75" Wide	31"
Seat depth	19.5"	19.75"	19.5"
Seat-to-back height	32"	32"	32"
Seat-to-floor height	21.5"	21"	22"
Overall tray width and length	25.5" W x 15" L	25.5" W x 15" L (adult model only)	10" W x 19" L
Carton dimensions L x W x H	40" L x 32" W x 29" H	37" L x 29" W x 29" H	40" L x 32" W x 29" H
Weight capacity	250 lb.	250 lb. Adult 400 lb. Wide	250 lb.
Product weight	73 lb.	86 lb.	89 lb.
Shipping weight	86 lb.	105 lb.	105 lb.
Limited warranty			
Chair	3 years	3 years	3 years
Upholstery	1 year	1 year	1 year
Gas cylinder	n/a	n/a	1 year

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2009 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Form No. 04-038
Rev. 02/09