




SWIVEL SHOWER CHAIR 2.0



Swivel Shower Chair 2.0: MHSCII

O.T. AWARD WINNING 2017

- Redesigned swivel and bearing mechanism for added strength and stability
- Swivels 360° and locks in 90° increments
- Swivel action reduces skin shear
- Removeable arm rests
- Height adjustable legs from 20"-25"
- Padded seat, back and arm rest
- Drainage holes for easy water exit
- Stainless steel pin is spring loaded and self locking
- 300 lbs weight capacity
- Weight - 10 lbs
- Rustproof, lightweight aluminum
- Tool free assembly
- Fits most bathtubs



PRODUCT DETAILS

SPECIFICATION

Rotation	360 ⁰
Seat	16" x 13"
Arm To Arm Width	20"
Height	20"-25" adjustable
Weight	10 lbs
Weight Capacity	300 lbs
Tool	Tool Free Assembly
Warranty	Limited One Year

Measurement of Swivel Shower Chair



MOBB Health Care Ltd.
116A Viceroy Rd. Concord,
ON L4K 2M1 Canada
Toll: 1-866-238-7190
Tel: 905-669-0348
www.mobbhhc.com

