

# Wave dimensions

| User dimensions (inches)                                                                                                                                                                                                                                                                                                                                                              | Z210 small       | Z220 medium                 | Z230 large       |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|-----------------------------|------------------|
| Height                                                                                                                                                                                                                                                                                                                                                                                | up to 46         | 40–56                       | 50–74            |
|  <b>Key user dimension: height</b><br>Select the appropriate bath chair by the user's overall height.<br>Choose the size that allows for growth.<br><br><b>Important:</b> Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load. |                  |                             |                  |
| Item dimensions (inches)                                                                                                                                                                                                                                                                                                                                                              | Z210 small       | Z220 medium                 | Z230 large       |
| Seat depth                                                                                                                                                                                                                                                                                                                                                                            | 10½              | 13                          | 17               |
| Seat width                                                                                                                                                                                                                                                                                                                                                                            | 15               | 15                          | 17½              |
| Back height                                                                                                                                                                                                                                                                                                                                                                           | 23               | 26½                         | 32               |
| Seat height                                                                                                                                                                                                                                                                                                                                                                           | 2–7              | 2–7                         | 2–7              |
| Seat height (incl. tub stand)                                                                                                                                                                                                                                                                                                                                                         | 12½ – 18½        | 12½ – 18½                   | 12½ – 18½        |
| Seat height (incl. shower stand)                                                                                                                                                                                                                                                                                                                                                      | 27½              | 27½                         | 27½              |
| Seat height (incl. tub transfer base)<br>with long leg inserts                                                                                                                                                                                                                                                                                                                        | 21–26<br>24½–29½ | 21–26<br>24½–29½            | 21–26<br>24½–29½ |
| Length of mesh area when flat (without calf rest)                                                                                                                                                                                                                                                                                                                                     | 35               | 41                          | 50½              |
| Length of mesh area when flat (with calf rest)                                                                                                                                                                                                                                                                                                                                        | 42½              | 49                          | 61               |
| Base width                                                                                                                                                                                                                                                                                                                                                                            | 18               | 18                          | 21               |
| Base length (legs down)<br>with legs flat                                                                                                                                                                                                                                                                                                                                             | 24<br>29½        | 24<br>29½                   | 24<br>29½        |
| Tub stand height (legs flat)                                                                                                                                                                                                                                                                                                                                                          | 8                | 8                           | 8                |
| Tub stand length (legs flat)                                                                                                                                                                                                                                                                                                                                                          | 39               | 39                          | 39               |
| Tub stand width (legs flat)                                                                                                                                                                                                                                                                                                                                                           | 17½              | 17½                         | 17½              |
| Shower stand height                                                                                                                                                                                                                                                                                                                                                                   | 22               | 22                          | 22               |
| Shower stand width                                                                                                                                                                                                                                                                                                                                                                    | 27               | 27                          | 27               |
| Shower stand length                                                                                                                                                                                                                                                                                                                                                                   | 34               | 34                          | 34               |
| Tub transfer base footprint                                                                                                                                                                                                                                                                                                                                                           | 28 x 30          | 28 x 30                     | 28 x 30          |
| Max. tub height for transfer base<br>with long leg inserts                                                                                                                                                                                                                                                                                                                            | 20<br>24         | 20<br>24                    | 20<br>24         |
| Maximum working load (lb)                                                                                                                                                                                                                                                                                                                                                             | 75               | 150                         | 250              |
| Back and seat angle adjustments                                                                                                                                                                                                                                                                                                                                                       |                  | 0°, 15°, 30°, 45°, 60°, 75° |                  |



All three sizes are available in blue and pink.

# 2021 order form

Use dimension chart to select appropriate size.

## Conversion kits

The tool-free conversion kit converts your small chair to medium as your child grows.

### Fabric without calf rest

Blue Z268  
Pink Z269

(Includes seat, back, cover and chest strap)



### Fabric with calf rest

Blue Z266  
Pink Z267

(Includes seat, back, cover and chest strap)



|                                                                                                      |                               | Z210 small                    | Z220 medium                   | Z230 large                    |
|------------------------------------------------------------------------------------------------------|-------------------------------|-------------------------------|-------------------------------|-------------------------------|
| * Chair                                                                                              |                               | Z211 <input type="checkbox"/> | Z221 <input type="checkbox"/> | Z231 <input type="checkbox"/> |
| * Fabric cover                                                                                       | CHOOSE                        |                               |                               |                               |
|                                                                                                      | Blue                          | Z244 <input type="checkbox"/> | Z246 <input type="checkbox"/> | Z248 <input type="checkbox"/> |
|                                                                                                      | Pink                          | Z245 <input type="checkbox"/> | Z247 <input type="checkbox"/> | Z249 <input type="checkbox"/> |
| * Chest strap                                                                                        | CHOOSE                        |                               |                               |                               |
|                                                                                                      | With lateral positioning      | Z212 <input type="checkbox"/> | Z222 <input type="checkbox"/> | Z232 <input type="checkbox"/> |
|                                                                                                      | Without lateral positioning   | Z213 <input type="checkbox"/> | Z223 <input type="checkbox"/> | Z233 <input type="checkbox"/> |
| Head blocks                                                                                          |                               | Z214 <input type="checkbox"/> | Z214 <input type="checkbox"/> | Z234 <input type="checkbox"/> |
| Leg straps                                                                                           |                               | Z215 <input type="checkbox"/> | Z225 <input type="checkbox"/> | Z235 <input type="checkbox"/> |
| Calf rest (Calf rest must be chosen at the time of the original purchase. It cannot be added later.) |                               | Z216 <input type="checkbox"/> | Z226 <input type="checkbox"/> | Z236 <input type="checkbox"/> |
| Tub stand                                                                                            |                               | Z217 <input type="checkbox"/> | Z217 <input type="checkbox"/> | Z217 <input type="checkbox"/> |
| Shower stand                                                                                         |                               | Z218 <input type="checkbox"/> | Z218 <input type="checkbox"/> | Z218 <input type="checkbox"/> |
| Tub transfer base (calf rest recommended)                                                            | CHOOSE                        |                               |                               |                               |
|                                                                                                      | Counterclockwise into tub     | Z227 <input type="checkbox"/> | Z227 <input type="checkbox"/> | Z227 <input type="checkbox"/> |
|                                                                                                      | Clockwise into tub            | Z228 <input type="checkbox"/> | Z228 <input type="checkbox"/> | Z228 <input type="checkbox"/> |
| Legs for tub transfer base                                                                           | CHOOSE                        |                               |                               |                               |
|                                                                                                      | Standard                      | Z224 <input type="checkbox"/> | Z224 <input type="checkbox"/> | Z224 <input type="checkbox"/> |
|                                                                                                      | Long (tub edge over 20" high) | Z238 <input type="checkbox"/> | Z238 <input type="checkbox"/> | Z238 <input type="checkbox"/> |
| Additional head blocks                                                                               |                               | Z214 <input type="checkbox"/> | Z214 <input type="checkbox"/> | Z234 <input type="checkbox"/> |
| Additional fabric cover                                                                              |                               | Z219 <input type="checkbox"/> | Z229 <input type="checkbox"/> | Z239 <input type="checkbox"/> |
| Additional chest strap                                                                               | CHOOSE                        |                               |                               |                               |
|                                                                                                      | With lateral positioning      | Z212 <input type="checkbox"/> | Z222 <input type="checkbox"/> | Z232 <input type="checkbox"/> |
|                                                                                                      | Without lateral positioning   | Z213 <input type="checkbox"/> | Z223 <input type="checkbox"/> | Z233 <input type="checkbox"/> |
| Total                                                                                                |                               | Z210                          | Z220                          | Z230                          |

\* Required