

Wave dimensions

User dimensions (inches)	Z210 small	Z220 medium	Z230 large
Height	up to 46	40–56	50–74
 Key user dimension: height Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth. Important: Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.			
Item dimensions (inches)	Z210 small	Z220 medium	Z230 large
Seat depth	10½	13	17
Seat width	15	15	17½
Back height	23	26½	32
Seat height	2–7	2–7	2–7
Seat height (incl. tub stand)	12½ – 18½	12½ – 18½	12½ – 18½
Seat height (incl. shower stand)	27½	27½	27½
Seat height (incl. tub transfer base) with long leg inserts	21–26 24½–29½	21–26 24½–29½	21–26 24½–29½
Length of mesh area when flat (without calf rest)	35	41	50½
Length of mesh area when flat (with calf rest)	42½	49	61
Base width	18	18	21
Base length (legs down) with legs flat	24 29½	24 29½	24 29½
Tub stand height (legs flat)	8	8	8
Tub stand length (legs flat)	39	39	39
Tub stand width (legs flat)	17½	17½	17½
Shower stand height	22	22	22
Shower stand width	27	27	27
Shower stand length	34	34	34
Tub transfer base footprint	28 x 30	28 x 30	28 x 30
Max. tub height for transfer base with long leg inserts	20 24	20 24	20 24
Maximum working load (lb)	75	150	250
Back and seat angle adjustments	0°, 15°, 30°, 45°, 60°, 75°		



All three sizes are available in blue and pink.

2021 order form

Use dimension chart to select appropriate size.

Conversion kits

The tool-free conversion kit converts your small chair to medium as your child grows.

Fabric without calf rest

Blue Z268
Pink Z269

(Includes seat, back, cover and chest strap)



Fabric with calf rest

Blue Z266
Pink Z267

(Includes seat, back, cover and chest strap)



		Z210 small	Z220 medium	Z230 large
* Chair		Z211 ☐	Z221 ☐	Z231 ☐
* Fabric cover	CHOOSE			
	Blue	Z244 ☐	Z246 ☐	Z248 ☐
	Pink	Z245 ☐	Z247 ☐	Z249 ☐
* Chest strap	CHOOSE			
	With lateral positioning	Z212 ☐	Z222 ☐	Z232 ☐
	Without lateral positioning	Z213 ☐	Z223 ☐	Z233 ☐
Head blocks		Z214 ☐	Z214 ☐	Z234 ☐
Leg straps		Z215 ☐	Z225 ☐	Z235 ☐
Calf rest (Calf rest must be chosen at the time of the original purchase. It cannot be added later.)		Z216 ☐	Z226 ☐	Z236 ☐
Tub stand		Z217 ☐	Z217 ☐	Z217 ☐
Shower stand		Z218 ☐	Z218 ☐	Z218 ☐
Tub transfer base (calf rest recommended)	CHOOSE			
	Counterclockwise into tub	Z227 ☐	Z227 ☐	Z227 ☐
	Clockwise into tub	Z228 ☐	Z228 ☐	Z228 ☐
Legs for tub transfer base	CHOOSE			
	Standard	Z224 ☐	Z224 ☐	Z224 ☐
	Long (tub edge over 20" high)	Z238 ☐	Z238 ☐	Z238 ☐
Additional head blocks		Z214 ☐	Z214 ☐	Z234 ☐
Additional fabric cover		Z219 ☐	Z229 ☐	Z239 ☐
Additional chest strap	CHOOSE			
	With lateral positioning	Z212 ☐	Z222 ☐	Z232 ☐
	Without lateral positioning	Z213 ☐	Z223 ☐	Z233 ☐
Total		Z210	Z220	Z230

* Required