



Molift Quick Raiser 2

Sit-to-stand hoist with electrically adjustable leg base

Molift Quick Raiser 2 with electrically adjustable leg base is compact, stable and easy to manoeuvre. The unique wheel position and near perfect weight distribution of the base provides a good turning radius and easy manoeuvrability even in narrow spaces.

The inclined straight column imitates the natural pattern of movement and fulfils the hoisting needs of users with some degree of weight bearing capacity. Molift Quick Raiser 2 has a safe working load of 160 kg.

The low base, 10 cm, fits under most beds and chairs while the electrically adjustable legs enables accessibility almost everywhere.



Electrically adjustable legs

A perfect alternative when the hoist is used together with chairs with no accessibility underneath, such as electric wheelchairs and recliners.



Adjustable knee pads

The knee pads are height adjustable and can be turned around for better fitting and to give alternative support for the legs.



Battery indicator

Battery indicator on the hand control showing when charging is needed.



Unique wheel positioning

The lightweight construction combined with the unique wheel positioning, provides exceptional maneuverability.

**Molift Quick Raiser 2****Item no.**

Chassis width

610 mm

V-shape arm
U-shape armM04200
M04400Visit etac.com to see the full range of accessories.**Max user weight (SWL)**

160 kg

Total weight

36 kg

DimensionsLeg height 100 mm
External width of wheelbase 1200 mm**Technical**Battery
Charge time
Hoisting speed
Rated performanceNiMH 14,4 V - 2,6 Ah
approx. 3 hrs
50 mm/sec (75 kg)
40 lifts (75 kg, 50 cm)max
160 kg
max
353 lbs

Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution.

For the latest news and continuously updated product information – visit: www.etac.com