

# Yeti S746

TRAILMASTER

## Mid-size scooter

The Eclipse Trailmaster series offer modern design, style, performance & features such as an LCD digital dashboard. This mid-size line of products are for individuals looking to access their local community for shopping, appointments and visiting with family and friends.

- > Trailmaster Yeti 4-wheel scooter
- > Digital LCD dashboard
- > Full light package
- > Rear mono shock suspension
- > 50AH batteries with 5A charger
- > Colours: red, blue or white



\$3,995.00 MSRP



### S746 Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb°
136KG/300LBS	52.5" x 24" x 48"	94KG/207LBS	1.6HP	10KPH/6.25MPH	2 x 50AH	10"	28KM/18M	150CM/59"	10°/16%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

[www.eclipsemedical.com](http://www.eclipsemedical.com)

#112019EN