


# Aeromesh® Bath Products

Users and caregivers alike love Bodypoint® bath products, made with our unique Aeromesh®, a silky, porous material that lets soap, water, dirt and air flow through freely while protecting skin. Padded for comfort, quick to dry, and compatible with all major shower chair brands.



Aeromesh® fabric is soft, lightweight, and porous

## Aeromesh® Bath Belts

Improve safety and stability in the bath without compromising comfort or hygiene. Can be worn on the chest, abdomen, hips or legs to support users while bathing or in the pool. Three styles, all length-adjustable with D-rings. Machine wash and dry.



### Two-piece pad with grommets

Center-release for quick & easy access. Stainless steel grommets for attachment.

**BB101M-1 Medium**  
61-94 cm (24-27") long

**BB101L-1 Large**  
79-109 cm (31-43") long



### One-piece pad with grommets

Side-release buckles for easy belt removal and cleaning. Mounting straps with stainless steel grommets.

**BB111L-1 Large**  
73.5-129.5 cm (29-51") long

**BB111XL-1 Extra Large**  
81-137 cm (32-54") long

Mounting hardware for grommet-style Bath Belts not included. Attach to chair frame with a stainless steel screw, Bodypoint Quick Release Knobs, or our non-rusting Band Clamp or Frame Clamp.



### One-piece pad, plastic slides

Wrap webbing around the bath chair frame and secure with plastic slides. Detach belt with side-release buckles.

**BB112L-1 Large**  
68.5 - 160 cm (27-63") long

**BB112XL-1 Extra Large**  
76-173 cm (30-68") long

## Shower Chair Calf Support

This clip-on calf support, designed specifically for use while bathing, helps stabilize the legs to reduce risk and increase comfort, especially for users with poor trunk control or in tilt position.



Clips fit on Ø22mm (7/8") and 25mm (1") leg rest tubing used by all major shower chair brands. Machine wash and dry.

**BB216-22MM** Fits 40.5 cm (16") wide chair  
**BB218-22MM** Fits 45.5 cm (18") wide chair  
**BB220-22MM** Fits 51 cm (20") wide chair

