

2000IPS Car Seat

FAX TO 905-787-8550 OR SUBMIT ONLINE



Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO #: _____ Ordered By: _____

Phone: _____ Fax: _____ Email: _____

Ship to Address: _____

Size	Standard Seat Depth	Standard Seat Depth w/ Abductor	Medium Seat Depth	Medium Seat Depth w/ Abductor	Large Seat Depth	Large Seat Depth w/ Abductor
Speedway Gray 	CSI-2000W-SG	CSI-2000WA-SG	CSI-2000WM-SG	CSI-2000WMA-SG	CSI-2000WL-SG	CSI-2000WLA-SG
SELECT TO ITEMIZE		CSI-2000WA-XX-IT	CSI-2000WM-XX-IT	CSI-2000WMA-XX-IT	CSI-2000WL-XX-IT	CSI-2000WLA-XX-IT
Seat Depth Extender			☒ CSI-K2025	☒ CSI-K2025	☒ CSI-K2026	☒ CSI-K2026
Abductor		☒ CSI-K2015		☒ CSI-K2018		☒ CSI-K2019
Weight Capacity	20–102 lb.	20–102 lb.	20–102 lb.	20–102 lb.	20–102 lb.	20–102 lb.
Max. Height Capacity	60"	60"	60"	60"	60"	60"
Harness Height	11.5", 14.5", 17.5", 20.5"					
Back Height	29"	29"	29"	29"	29"	29"
Seat Depth	10"	10"	12.5"	12.5"	15"	15"
Inside Seat Width, Front	15.5"	15.5"	15.5"	15.5"	15.5"	15.5"
Inside Seat Width, Rear	13"	13"	13"	13"	13"	13"
Overall Width	19"	19"	19"	19"	19"	19"
Overall Depth (w/ seat rails)	13.5"	16.5"	16"	19"	18.5"	21.5"
Seat-to-Back Angle	80°	80°	80°	80°	80°	80°
Product Weight	16 lb.	18 lb.	19 lb.	22 lb.	20 lb.	25 lb.

ACCESSORIES (for all IPS Car Seats)

Block-Style Trunk Supports, Black Mesh CSI-A2008-BK

Lateral Pads CSI-A2010

Bolt-in Hardware Kit with Tether Strap CSI-A2031

Extensor Thrust Wedge, 2" CSI-A2050

Buckle Guard ZZ-A2075

Retainer Clip Guard for car seats CSZ-A2085